

С-Реактивный белок, (CRP-Uhs) Высококчувствительный.

ИтипоTD (Иммунотурбодиметрический с латексными частицами)



Параметры для ввода в программу анализатора iLab-Taurus

Test Name: NN **CRPh** Test Code: IL_ **CRPh-d** **Measure** **Sample Reagent** **Ranges Limits** **Calibration**

Sample Type
Serum

Reaction Cycle
☒ Standard ☐ Extended

Reporting Unit
☒ Мг/л* ☐ User Define
 Decimal Points: 2
 Conv.Factor: 0.000

Methodology
Type
☐ End Point ☒ Rate
Measuring Point
 19 - ☒ 29
Photometric
☒ 1 Wavelength ☐ 2 Wavelength
Primary: 510 **Secondary**:
Correction Constant
Slope: 1.000 **Intercept**: 0.000

Test Name: NN **CRPh** Test Code: IL_ **CRPh-d** **Measure** **Sample Reagent** **Ranges Limits** **Calibration**

Sample Volume

	Sample Volume	Dilution Sample Vol.	Dilution Diluent Vol.
1	10.0	0.0	0.0
2	2.0	0.0	0.0
3	12.0	0.0	0.0
4	0.0	0.0	0.0

Diluent:
Diluent Warning Limit: 0 Tests

Reagent Volume

R1: **CRPh-d1** R2: **CRPh-d2**

	Volume	Diluent Vol.	stirring	Warning Limit(tests)	Stability(days)
R1	100	0.0	<input checked="" type="checkbox"/>	20	No Control
R2	100	0.0	<input checked="" type="checkbox"/>	20	No Control

Sampling Condition

Condition No.	1	2	3	4
First Run	<input checked="" type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
Samp.Vol.Reduction	<input type="radio"/>	<input checked="" type="radio"/>	<input type="radio"/>	<input type="radio"/>
Below N-Range	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
Above N-Range	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
Panic L	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
Panic H	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
User Range L	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
<input checked="" type="checkbox"/> User Range H	<input type="radio"/>	<input checked="" type="radio"/>	<input type="radio"/>	<input type="radio"/>
Noise	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
Prozone	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
<input checked="" type="checkbox"/> HIGH!	<input type="radio"/>	<input checked="" type="radio"/>	<input type="radio"/>	<input type="radio"/>
<input checked="" type="checkbox"/> ABS!	<input type="radio"/>	<input checked="" type="radio"/>	<input type="radio"/>	<input type="radio"/>

C-Реактивный белок, (CRP-Uhs) Высококчувствительный.

ИтипоTD (Иммунотурбодиметрический с латексными частицами)



Ranges		Limits																																																																																																		
<table border="1"> <thead> <tr> <th colspan="3">Normal Range</th> </tr> <tr> <th></th> <th>Lower</th> <th>Upper</th> </tr> </thead> <tbody> <tr> <td>Male</td> <td>0.0</td> <td>5.0</td> </tr> <tr> <td>Female</td> <td>0.0</td> <td>5.0</td> </tr> <tr> <td>Other</td> <td>0.0</td> <td>5.0</td> </tr> </tbody> </table>		Normal Range				Lower	Upper	Male	0.0	5.0	Female	0.0	5.0	Other	0.0	5.0	<table border="1"> <thead> <tr> <th colspan="2">Reaction Slope</th> <th colspan="2">Absorbance Limit</th> </tr> </thead> <tbody> <tr> <td><input type="radio"/> Negative</td> <td><input checked="" type="radio"/> Positive</td> <td><input checked="" type="radio"/> Above</td> <td><input type="radio"/> Below</td> </tr> <tr> <td colspan="2"> <table border="1"> <thead> <tr> <th colspan="2">Non-Linear Limit</th> </tr> <tr> <td>100</td> <td>%</td> </tr> </thead> <tbody> <tr> <td><input type="text"/></td> <td><input type="text"/></td> </tr> </tbody> </table> </td> <td colspan="2"> <table border="1"> <thead> <tr> <th colspan="2">Prozone Limit</th> </tr> <tr> <td><input type="radio"/> Above</td> <td><input checked="" type="radio"/> Below</td> </tr> </thead> <tbody> <tr> <td> <table border="1"> <thead> <tr> <th colspan="2">Limit</th> </tr> <tr> <td>0.0</td> <td><input type="text"/></td> </tr> </thead> <tbody> <tr> <td><input type="text"/></td> <td><input type="text"/></td> </tr> </tbody> </table> </td> <td> <table border="1"> <thead> <tr> <th colspan="2">Equation</th> </tr> <tr> <td>none</td> <td><input type="text"/></td> </tr> </thead> <tbody> <tr> <td><input type="text"/></td> <td><input type="text"/></td> </tr> </tbody> </table> </td> </tr> <tr> <td colspan="2"> <table border="1"> <thead> <tr> <th colspan="2">Judge Point</th> </tr> <tr> <td>0</td> <td><input type="text"/></td> </tr> </thead> <tbody> <tr> <td><input type="text"/></td> <td><input type="text"/></td> </tr> </tbody> </table> </td> <td colspan="2"></td> </tr> </tbody> </table> </td> <td colspan="2"> <table border="1"> <thead> <tr> <th colspan="2">User Range</th> </tr> <tr> <th>Lower</th> <th>Upper</th> </tr> </thead> <tbody> <tr> <td>0.01</td> <td>35</td> </tr> </tbody> </table> </td> <td colspan="2"> <table border="1"> <thead> <tr> <th colspan="2">ValidRange</th> </tr> <tr> <th>Lower</th> <th>Upper</th> </tr> </thead> <tbody> <tr> <td>0.01</td> <td>35</td> </tr> </tbody> </table> </td> </tr> <tr> <td colspan="2"> <table border="1"> <thead> <tr> <th colspan="2">Qualitative</th> </tr> <tr> <td><input type="radio"/> On</td> <td><input checked="" type="radio"/> Off</td> </tr> </thead> <tbody> <tr> <td><input type="text"/></td> <td><input type="text"/></td> </tr> </tbody> </table> </td> <td colspan="2"> <table border="1"> <thead> <tr> <th colspan="2">Serum Index Limits</th> </tr> <tr> <th></th> <th>Limit</th> </tr> </thead> <tbody> <tr> <td>Hemolysis</td> <td>0.0</td> </tr> <tr> <td>Icterus</td> <td>0.0</td> </tr> <tr> <td>Lipemia</td> <td>0.0</td> </tr> </tbody> </table> </td> </tr> </tbody> </table>		Reaction Slope		Absorbance Limit		<input type="radio"/> Negative	<input checked="" type="radio"/> Positive	<input checked="" type="radio"/> Above	<input type="radio"/> Below	<table border="1"> <thead> <tr> <th colspan="2">Non-Linear Limit</th> </tr> <tr> <td>100</td> <td>%</td> </tr> </thead> <tbody> <tr> <td><input type="text"/></td> <td><input type="text"/></td> </tr> </tbody> </table>		Non-Linear Limit		100	%	<input type="text"/>	<input type="text"/>	<table border="1"> <thead> <tr> <th colspan="2">Prozone Limit</th> </tr> <tr> <td><input type="radio"/> Above</td> <td><input checked="" type="radio"/> Below</td> </tr> </thead> <tbody> <tr> <td> <table border="1"> <thead> <tr> <th colspan="2">Limit</th> </tr> <tr> <td>0.0</td> <td><input type="text"/></td> </tr> </thead> <tbody> <tr> <td><input type="text"/></td> <td><input type="text"/></td> </tr> </tbody> </table> </td> <td> <table border="1"> <thead> <tr> <th colspan="2">Equation</th> </tr> <tr> <td>none</td> <td><input type="text"/></td> </tr> </thead> <tbody> <tr> <td><input type="text"/></td> <td><input type="text"/></td> </tr> </tbody> </table> </td> </tr> <tr> <td colspan="2"> <table border="1"> <thead> <tr> <th colspan="2">Judge Point</th> </tr> <tr> <td>0</td> <td><input type="text"/></td> </tr> </thead> <tbody> <tr> <td><input type="text"/></td> <td><input type="text"/></td> </tr> </tbody> </table> </td> <td colspan="2"></td> </tr> </tbody> </table>		Prozone Limit		<input type="radio"/> Above	<input checked="" type="radio"/> Below	<table border="1"> <thead> <tr> <th colspan="2">Limit</th> </tr> <tr> <td>0.0</td> <td><input type="text"/></td> </tr> </thead> <tbody> <tr> <td><input type="text"/></td> <td><input type="text"/></td> </tr> </tbody> </table>	Limit		0.0	<input type="text"/>	<input type="text"/>	<input type="text"/>	<table border="1"> <thead> <tr> <th colspan="2">Equation</th> </tr> <tr> <td>none</td> <td><input type="text"/></td> </tr> </thead> <tbody> <tr> <td><input type="text"/></td> <td><input type="text"/></td> </tr> </tbody> </table>	Equation		none	<input type="text"/>	<input type="text"/>	<input type="text"/>	<table border="1"> <thead> <tr> <th colspan="2">Judge Point</th> </tr> <tr> <td>0</td> <td><input type="text"/></td> </tr> </thead> <tbody> <tr> <td><input type="text"/></td> <td><input type="text"/></td> </tr> </tbody> </table>		Judge Point		0	<input type="text"/>	<input type="text"/>	<input type="text"/>			<table border="1"> <thead> <tr> <th colspan="2">User Range</th> </tr> <tr> <th>Lower</th> <th>Upper</th> </tr> </thead> <tbody> <tr> <td>0.01</td> <td>35</td> </tr> </tbody> </table>		User Range		Lower	Upper	0.01	35	<table border="1"> <thead> <tr> <th colspan="2">ValidRange</th> </tr> <tr> <th>Lower</th> <th>Upper</th> </tr> </thead> <tbody> <tr> <td>0.01</td> <td>35</td> </tr> </tbody> </table>		ValidRange		Lower	Upper	0.01	35	<table border="1"> <thead> <tr> <th colspan="2">Qualitative</th> </tr> <tr> <td><input type="radio"/> On</td> <td><input checked="" type="radio"/> Off</td> </tr> </thead> <tbody> <tr> <td><input type="text"/></td> <td><input type="text"/></td> </tr> </tbody> </table>		Qualitative		<input type="radio"/> On	<input checked="" type="radio"/> Off	<input type="text"/>	<input type="text"/>	<table border="1"> <thead> <tr> <th colspan="2">Serum Index Limits</th> </tr> <tr> <th></th> <th>Limit</th> </tr> </thead> <tbody> <tr> <td>Hemolysis</td> <td>0.0</td> </tr> <tr> <td>Icterus</td> <td>0.0</td> </tr> <tr> <td>Lipemia</td> <td>0.0</td> </tr> </tbody> </table>		Serum Index Limits			Limit	Hemolysis	0.0	Icterus	0.0	Lipemia	0.0
Normal Range																																																																																																				
	Lower	Upper																																																																																																		
Male	0.0	5.0																																																																																																		
Female	0.0	5.0																																																																																																		
Other	0.0	5.0																																																																																																		
Reaction Slope		Absorbance Limit																																																																																																		
<input type="radio"/> Negative	<input checked="" type="radio"/> Positive	<input checked="" type="radio"/> Above	<input type="radio"/> Below																																																																																																	
<table border="1"> <thead> <tr> <th colspan="2">Non-Linear Limit</th> </tr> <tr> <td>100</td> <td>%</td> </tr> </thead> <tbody> <tr> <td><input type="text"/></td> <td><input type="text"/></td> </tr> </tbody> </table>		Non-Linear Limit		100	%	<input type="text"/>	<input type="text"/>	<table border="1"> <thead> <tr> <th colspan="2">Prozone Limit</th> </tr> <tr> <td><input type="radio"/> Above</td> <td><input checked="" type="radio"/> Below</td> </tr> </thead> <tbody> <tr> <td> <table border="1"> <thead> <tr> <th colspan="2">Limit</th> </tr> <tr> <td>0.0</td> <td><input type="text"/></td> </tr> </thead> <tbody> <tr> <td><input type="text"/></td> <td><input type="text"/></td> </tr> </tbody> </table> </td> <td> <table border="1"> <thead> <tr> <th colspan="2">Equation</th> </tr> <tr> <td>none</td> <td><input type="text"/></td> </tr> </thead> <tbody> <tr> <td><input type="text"/></td> <td><input type="text"/></td> </tr> </tbody> </table> </td> </tr> <tr> <td colspan="2"> <table border="1"> <thead> <tr> <th colspan="2">Judge Point</th> </tr> <tr> <td>0</td> <td><input type="text"/></td> </tr> </thead> <tbody> <tr> <td><input type="text"/></td> <td><input type="text"/></td> </tr> </tbody> </table> </td> <td colspan="2"></td> </tr> </tbody> </table>		Prozone Limit		<input type="radio"/> Above	<input checked="" type="radio"/> Below	<table border="1"> <thead> <tr> <th colspan="2">Limit</th> </tr> <tr> <td>0.0</td> <td><input type="text"/></td> </tr> </thead> <tbody> <tr> <td><input type="text"/></td> <td><input type="text"/></td> </tr> </tbody> </table>	Limit		0.0	<input type="text"/>	<input type="text"/>	<input type="text"/>	<table border="1"> <thead> <tr> <th colspan="2">Equation</th> </tr> <tr> <td>none</td> <td><input type="text"/></td> </tr> </thead> <tbody> <tr> <td><input type="text"/></td> <td><input type="text"/></td> </tr> </tbody> </table>	Equation		none	<input type="text"/>	<input type="text"/>	<input type="text"/>	<table border="1"> <thead> <tr> <th colspan="2">Judge Point</th> </tr> <tr> <td>0</td> <td><input type="text"/></td> </tr> </thead> <tbody> <tr> <td><input type="text"/></td> <td><input type="text"/></td> </tr> </tbody> </table>		Judge Point		0	<input type="text"/>	<input type="text"/>	<input type="text"/>			<table border="1"> <thead> <tr> <th colspan="2">User Range</th> </tr> <tr> <th>Lower</th> <th>Upper</th> </tr> </thead> <tbody> <tr> <td>0.01</td> <td>35</td> </tr> </tbody> </table>		User Range		Lower	Upper	0.01	35	<table border="1"> <thead> <tr> <th colspan="2">ValidRange</th> </tr> <tr> <th>Lower</th> <th>Upper</th> </tr> </thead> <tbody> <tr> <td>0.01</td> <td>35</td> </tr> </tbody> </table>		ValidRange		Lower	Upper	0.01	35																																															
Non-Linear Limit																																																																																																				
100	%																																																																																																			
<input type="text"/>	<input type="text"/>																																																																																																			
Prozone Limit																																																																																																				
<input type="radio"/> Above	<input checked="" type="radio"/> Below																																																																																																			
<table border="1"> <thead> <tr> <th colspan="2">Limit</th> </tr> <tr> <td>0.0</td> <td><input type="text"/></td> </tr> </thead> <tbody> <tr> <td><input type="text"/></td> <td><input type="text"/></td> </tr> </tbody> </table>	Limit		0.0	<input type="text"/>	<input type="text"/>	<input type="text"/>	<table border="1"> <thead> <tr> <th colspan="2">Equation</th> </tr> <tr> <td>none</td> <td><input type="text"/></td> </tr> </thead> <tbody> <tr> <td><input type="text"/></td> <td><input type="text"/></td> </tr> </tbody> </table>	Equation		none	<input type="text"/>	<input type="text"/>	<input type="text"/>																																																																																							
Limit																																																																																																				
0.0	<input type="text"/>																																																																																																			
<input type="text"/>	<input type="text"/>																																																																																																			
Equation																																																																																																				
none	<input type="text"/>																																																																																																			
<input type="text"/>	<input type="text"/>																																																																																																			
<table border="1"> <thead> <tr> <th colspan="2">Judge Point</th> </tr> <tr> <td>0</td> <td><input type="text"/></td> </tr> </thead> <tbody> <tr> <td><input type="text"/></td> <td><input type="text"/></td> </tr> </tbody> </table>		Judge Point		0	<input type="text"/>	<input type="text"/>	<input type="text"/>																																																																																													
Judge Point																																																																																																				
0	<input type="text"/>																																																																																																			
<input type="text"/>	<input type="text"/>																																																																																																			
User Range																																																																																																				
Lower	Upper																																																																																																			
0.01	35																																																																																																			
ValidRange																																																																																																				
Lower	Upper																																																																																																			
0.01	35																																																																																																			
<table border="1"> <thead> <tr> <th colspan="2">Qualitative</th> </tr> <tr> <td><input type="radio"/> On</td> <td><input checked="" type="radio"/> Off</td> </tr> </thead> <tbody> <tr> <td><input type="text"/></td> <td><input type="text"/></td> </tr> </tbody> </table>		Qualitative		<input type="radio"/> On	<input checked="" type="radio"/> Off	<input type="text"/>	<input type="text"/>	<table border="1"> <thead> <tr> <th colspan="2">Serum Index Limits</th> </tr> <tr> <th></th> <th>Limit</th> </tr> </thead> <tbody> <tr> <td>Hemolysis</td> <td>0.0</td> </tr> <tr> <td>Icterus</td> <td>0.0</td> </tr> <tr> <td>Lipemia</td> <td>0.0</td> </tr> </tbody> </table>		Serum Index Limits			Limit	Hemolysis	0.0	Icterus	0.0	Lipemia	0.0																																																																																	
Qualitative																																																																																																				
<input type="radio"/> On	<input checked="" type="radio"/> Off																																																																																																			
<input type="text"/>	<input type="text"/>																																																																																																			
Serum Index Limits																																																																																																				
	Limit																																																																																																			
Hemolysis	0.0																																																																																																			
Icterus	0.0																																																																																																			
Lipemia	0.0																																																																																																			

Test Name		Test Code		Measure		Sample Reagent		Ranges Limits		Calibration																																																																																	
NN	CRPh	IL	CRPh-d																																																																																								
Calibration				Condition				Calibration/R-Blank Limit																																																																																			
<table border="1"> <thead> <tr> <th colspan="2">Method</th> <th colspan="2">Curve Type</th> </tr> </thead> <tbody> <tr> <td>M-point</td> <td><input type="text"/></td> <td>5P Log</td> <td><input type="text"/></td> </tr> </tbody> </table>				Method		Curve Type		M-point	<input type="text"/>	5P Log	<input type="text"/>	<table border="1"> <thead> <tr> <th colspan="2">Repeats</th> <th colspan="2">Stability</th> </tr> </thead> <tbody> <tr> <td>2</td> <td><input type="text"/></td> <td></td> <td>days</td> </tr> </tbody> </table>				Repeats		Stability		2	<input type="text"/>		days	<table border="1"> <thead> <tr> <th colspan="2">No</th> </tr> </thead> <tbody> <tr> <td>6</td> <td><input type="text"/></td> </tr> </tbody> </table>				No		6	<input type="text"/>	<table border="1"> <thead> <tr> <th colspan="2">Calibrator</th> </tr> <tr> <th></th> <th>Conc</th> </tr> </thead> <tbody> <tr> <td>1</td> <td>NaCl-0.9</td> </tr> <tr> <td>2</td> <td>TC-CRPh-1</td> </tr> <tr> <td>3</td> <td>TC-CRPh-2</td> </tr> <tr> <td>4</td> <td>TC-CRPh-3</td> </tr> <tr> <td>5</td> <td>TC-CRPh-4</td> </tr> <tr> <td>6</td> <td>TC-CRPh-5</td> </tr> </tbody> </table>				Calibrator			Conc	1	NaCl-0.9	2	TC-CRPh-1	3	TC-CRPh-2	4	TC-CRPh-3	5	TC-CRPh-4	6	TC-CRPh-5	<table border="1"> <thead> <tr> <th colspan="2">R-Blank Limit</th> </tr> </thead> <tbody> <tr> <td>2500.0</td> <td>mAbs</td> </tr> </tbody> </table>				R-Blank Limit		2500.0	mAbs	<table border="1"> <thead> <tr> <th colspan="2">Cal Reps Range</th> </tr> </thead> <tbody> <tr> <td>20.0</td> <td>%</td> </tr> </tbody> </table>				Cal Reps Range		20.0	%	<table border="1"> <thead> <tr> <th colspan="2">Min Cal Reps</th> </tr> </thead> <tbody> <tr> <td>100.0</td> <td>mAbs</td> </tr> </tbody> </table>				Min Cal Reps		100.0	mAbs	<table border="1"> <thead> <tr> <th colspan="2">Factor Change</th> </tr> </thead> <tbody> <tr> <td>35</td> <td>%</td> </tr> </tbody> </table>				Factor Change		35	%	<table border="1"> <thead> <tr> <th colspan="2">M-Point Curve Fit</th> </tr> </thead> <tbody> <tr> <td>0.0</td> <td>%</td> </tr> </tbody> </table>				M-Point Curve Fit		0.0	%
Method		Curve Type																																																																																									
M-point	<input type="text"/>	5P Log	<input type="text"/>																																																																																								
Repeats		Stability																																																																																									
2	<input type="text"/>		days																																																																																								
No																																																																																											
6	<input type="text"/>																																																																																										
Calibrator																																																																																											
	Conc																																																																																										
1	NaCl-0.9																																																																																										
2	TC-CRPh-1																																																																																										
3	TC-CRPh-2																																																																																										
4	TC-CRPh-3																																																																																										
5	TC-CRPh-4																																																																																										
6	TC-CRPh-5																																																																																										
R-Blank Limit																																																																																											
2500.0	mAbs																																																																																										
Cal Reps Range																																																																																											
20.0	%																																																																																										
Min Cal Reps																																																																																											
100.0	mAbs																																																																																										
Factor Change																																																																																											
35	%																																																																																										
M-Point Curve Fit																																																																																											
0.0	%																																																																																										

С-Реактивный белок, (CRP-Uhs) Высокочувствительный.

ИтипоTD (Иммунотурбодиметрический с латексными частицами)



* -Рекомендуемый параметр.

** - Вводится из паспорта калибратору.

Калибровка многоточечная (6 точек, 0 + 5 ур.) Spline или 5P Log, калибратор TruCal CRP-hs.

В качестве нулевой точки используйте физраствор.

Контроль TruLab CRP-hs уровень 1 и 2.