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Test Name: NN **GGT** Test Code: IL_ **GGT-d** Measure Sample Reagent Ranges Limits Calibration

Sample Type: **Serum**

Reaction Cycle: ☒ Standard ☐ Extended

Reporting Unit: ☒ Ед/л * ☐ User Define
 Decimal Points: 1
 Conv. Factor: 0.000

Methodology:
 Type: ☐ End Point ☒ Rate
 Measuring Point: 19 - ☒ 29
 Photometric: ☐ 1 Wavelength ☒ 2 Wavelength
 Primary: 405 Secondary: 546
 Correction Constant: Slope: 1.000 Intercept: 0.000

Test Name: NN **GGT** Test Code: IL_ **GGT-d** Measure Sample Reagent Ranges Limits Calibration

Sample Volume:

	Sample Volume	Dilution Sample Vol.	Dilution Diluent Vol.
1	14.0	0.0	0.0
2	2.0	0.0	0.0
3	25.0	0.0	0.0
4	10.0	5.0	95.0

Diluent:
 Diluent Warning Limit: 0 Tests

Reagent Volume:
 R1: **GGT-d1** R2: **GGT-d2**

	Volume	Diluent Vol.	stirring	Warning Limit(tests)	Stability(days)
R1	140	0.0	<input checked="" type="checkbox"/>	20	No Control
R2	35	0.0	<input checked="" type="checkbox"/>	20	No Control

Sampling Condition:

Condition No.	1	2	3	4
First Run	<input checked="" type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
Samp.Vol.Reduction	<input type="radio"/>	<input checked="" type="radio"/>	<input type="radio"/>	<input type="radio"/>
Below N-Range	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
Above N-Range	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
Panic L	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
Panic H	<input type="radio"/>	<input checked="" type="radio"/>	<input type="radio"/>	<input type="radio"/>
User Range L	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
<input checked="" type="checkbox"/> User Range H	<input type="radio"/>	<input checked="" type="radio"/>	<input type="radio"/>	<input type="radio"/>
Noise	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
Prozone	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
<input checked="" type="checkbox"/> HIGH!	<input type="radio"/>	<input checked="" type="radio"/>	<input type="radio"/>	<input type="radio"/>
<input checked="" type="checkbox"/> ABS!	<input type="radio"/>	<input checked="" type="radio"/>	<input type="radio"/>	<input type="radio"/>

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Ranges			Limits													
Normal Range <table border="1"> <thead> <tr> <th></th> <th>Lower</th> <th>Upper</th> </tr> </thead> <tbody> <tr> <td>Male</td> <td>0</td> <td>49</td> </tr> <tr> <td>Female</td> <td>0</td> <td>32</td> </tr> <tr> <td>Other</td> <td>0</td> <td>49</td> </tr> </tbody> </table>				Lower	Upper	Male	0	49	Female	0	32	Other	0	49	Reaction Slope <input type="radio"/> Negative <input checked="" type="radio"/> Positive	
	Lower	Upper														
Male	0	49														
Female	0	32														
Other	0	49														
User Range <table border="1"> <thead> <tr> <th>Lower</th> <th>Upper</th> </tr> </thead> <tbody> <tr> <td>2</td> <td>1200</td> </tr> </tbody> </table>			Lower	Upper	2	1200	Absorbance Limit <input checked="" type="radio"/> Above <input type="radio"/> Below 2500.0 mAbs									
Lower	Upper															
2	1200															
Valid Range <table border="1"> <thead> <tr> <th>Lower</th> <th>Upper</th> </tr> </thead> <tbody> <tr> <td>2</td> <td>1200</td> </tr> </tbody> </table>			Lower	Upper	2	1200	Non-Linear Limit 25 %									
Lower	Upper															
2	1200															
Qualitative <input type="radio"/> On <input checked="" type="radio"/> Off			Prozone Limit <input type="radio"/> Above <input checked="" type="radio"/> Below Limit: 0.0 Equation: none Judge Point: 0													
			Serum Index Limits <table border="1"> <thead> <tr> <th></th> <th>Limit</th> </tr> </thead> <tbody> <tr> <td>Hemolysis</td> <td>0.0</td> </tr> <tr> <td>Icterus</td> <td>0.0</td> </tr> <tr> <td>Lipemia</td> <td>0.0</td> </tr> </tbody> </table>			Limit	Hemolysis	0.0	Icterus	0.0	Lipemia	0.0				
	Limit															
Hemolysis	0.0															
Icterus	0.0															
Lipemia	0.0															

Test Name		Test Code		Measure	Sample Reagent	Ranges Limits	Calibration																		
NN	GGT	IL_	GGT-d																						
Calibration		Condition			Calibration/R-Blank Limit																				
Method 1-point Curve Type Linear Repeats 2 No Stability days		<input checked="" type="checkbox"/> Reagent Blank <input type="checkbox"/> Auto Reagent Blank by Bottle <input type="checkbox"/> Ask for calibration when reagent lot changes			R-Blank Limit 2500.0 mAbs Cal Reps Range 20.0 % Min Cal Reps 50.0 mAbs Factor Change 20 % M-Point Curve Fit 0.0 %																				
Calibrator <table border="1"> <thead> <tr> <th></th> <th>Calibrator</th> <th>Conc</th> </tr> </thead> <tbody> <tr> <td>1</td> <td>TruCalU</td> <td>**</td> </tr> <tr> <td>2</td> <td></td> <td></td> </tr> <tr> <td>3</td> <td></td> <td></td> </tr> <tr> <td>4</td> <td></td> <td></td> </tr> <tr> <td>5</td> <td></td> <td></td> </tr> </tbody> </table>			Calibrator	Conc	1	TruCalU	**	2			3			4			5								
	Calibrator	Conc																							
1	TruCalU	**																							
2																									
3																									
4																									
5																									

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