

Параметры для ввода в программу анализатора iLab-Taurus

Test Name: NN **Lipase** Test Code: IL_ **Lps -d** **Measure** **Sample Reagent** **Ranges Limits** **Calibration**

Sample Type
Serum

Reaction Cycle
☒ Standard ☐ Extended

Reporting Unit
☒ Ед/л * ☐ User Define
Decimal Points: 0
Conv. Factor: 0.000

Methodology
Type
☐ End Point ☒ Rate
Measuring Point
22 - ☒ 31
Photometric
☐ 1 Wavelength ☒ 2 Wavelength
Primary 570 **Secondary** 700
Correction Constant
Slope 1.000 **Intercept** 0.000

Test Name: NN **Lipase** Test Code: IL_ **Lps -d** **Measure** **Sample Reagent** **Ranges Limits** **Calibration**

Sample Volume

	Sample	Dilution	
	Volume	Sample Vol.	Diluent Vol.
1	3.25	0.0	0.0
2	3.0	10.0	90.0
3	6.0	0.0	0.0
4	2.0	5.0	95.0

Diluent
Diluent Warning Limit 0 Tests

Reagent Volume

R1: aAmy-d1 R2: aAmy-d2

	Volume	Diluent Vol.	stirring	Warning Limit(tests)	Stability(days)
R1	160	0.0	<input checked="" type="checkbox"/>	20	No Control
R2	40	0.0	<input checked="" type="checkbox"/>	20	No Control

Sampling Condition

	Condition No.	1	2	3	4
<input type="checkbox"/>	First Run	<input checked="" type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
<input type="checkbox"/>	Samp.Vol.Reduction	<input type="radio"/>	<input checked="" type="radio"/>	<input type="radio"/>	<input type="radio"/>
<input type="checkbox"/>	Below N-Range	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
<input type="checkbox"/>	Above N-Range	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
<input type="checkbox"/>	Panic L	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
<input type="checkbox"/>	Panic H	<input type="radio"/>	<input checked="" type="radio"/>	<input type="radio"/>	<input type="radio"/>
<input type="checkbox"/>	User Range L	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
<input checked="" type="checkbox"/>	User Range H	<input type="radio"/>	<input checked="" type="radio"/>	<input type="radio"/>	<input type="radio"/>
<input type="checkbox"/>	Noise	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
<input type="checkbox"/>	Prozone	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
<input checked="" type="checkbox"/>	HIGH!	<input type="radio"/>	<input checked="" type="radio"/>	<input type="radio"/>	<input type="radio"/>
<input checked="" type="checkbox"/>	ABS!	<input type="radio"/>	<input checked="" type="radio"/>	<input type="radio"/>	<input type="radio"/>

Ranges			Limits																																																																																																			
<table border="1"> <thead> <tr> <th colspan="3">Normal Range</th> </tr> <tr> <th></th> <th>Lower</th> <th>Upper</th> </tr> </thead> <tbody> <tr> <td>Male</td> <td>0</td> <td>60</td> </tr> <tr> <td>Female</td> <td>0</td> <td>60</td> </tr> <tr> <td>Other</td> <td>0</td> <td>60</td> </tr> </tbody> </table>			Normal Range				Lower	Upper	Male	0	60	Female	0	60	Other	0	60	<table border="1"> <thead> <tr> <th colspan="2">Reaction Slope</th> <th colspan="2">Absorbance Limit</th> </tr> </thead> <tbody> <tr> <td><input type="radio"/> Negative</td> <td><input checked="" type="radio"/> Positive</td> <td><input checked="" type="radio"/> Above</td> <td><input type="radio"/> Below</td> </tr> <tr> <td colspan="2"> <table border="1"> <thead> <tr> <th colspan="2">Non-Linear Limit</th> </tr> <tr> <th>Value</th> <th>%</th> </tr> </thead> <tbody> <tr> <td>25</td> <td></td> </tr> </tbody> </table> </td> <td colspan="2"> <table border="1"> <thead> <tr> <th colspan="2">Prozone Limit</th> </tr> <tr> <th>Value</th> <th>%</th> </tr> </thead> <tbody> <tr> <td>2500.0</td> <td></td> </tr> </tbody> </table> </td> </tr> <tr> <td colspan="2"> <table border="1"> <thead> <tr> <th colspan="2">User Range</th> </tr> <tr> <th>Lower</th> <th>Upper</th> </tr> </thead> <tbody> <tr> <td>3</td> <td>300</td> </tr> </tbody> </table> </td> <td colspan="2"> <table border="1"> <thead> <tr> <th colspan="2">ValidRange</th> </tr> <tr> <th>Lower</th> <th>Upper</th> </tr> </thead> <tbody> <tr> <td>3</td> <td>300</td> </tr> </tbody> </table> </td> <td colspan="2"> <table border="1"> <thead> <tr> <th colspan="2">Serum Index Limits</th> </tr> <tr> <th>Index</th> <th>Limit</th> </tr> </thead> <tbody> <tr> <td>Hemolysis</td> <td>0.0</td> </tr> <tr> <td>Icterus</td> <td>0.0</td> </tr> <tr> <td>Lipemia</td> <td>0.0</td> </tr> </tbody> </table> </td> </tr> <tr> <td colspan="3"> <table border="1"> <thead> <tr> <th colspan="2">Qualitative</th> </tr> <tr> <th>On</th> <th>Off</th> </tr> </thead> <tbody> <tr> <td><input type="radio"/></td> <td><input checked="" type="radio"/></td> </tr> </tbody> </table> </td> <td colspan="3"> <table border="1"> <thead> <tr> <th colspan="2">Limit</th> <th colspan="2">Equation</th> </tr> </thead> <tbody> <tr> <td>0.0</td> <td></td> <td>none</td> <td></td> </tr> <tr> <td colspan="2"> <table border="1"> <thead> <tr> <th colspan="2">Judge Point</th> </tr> <tr> <th>Value</th> <th>%</th> </tr> </thead> <tbody> <tr> <td>0</td> <td></td> </tr> </tbody> </table> </td> <td colspan="2"></td> </tr> </tbody> </table> </td> </tr> </tbody> </table>			Reaction Slope		Absorbance Limit		<input type="radio"/> Negative	<input checked="" type="radio"/> Positive	<input checked="" type="radio"/> Above	<input type="radio"/> Below	<table border="1"> <thead> <tr> <th colspan="2">Non-Linear Limit</th> </tr> <tr> <th>Value</th> <th>%</th> </tr> </thead> <tbody> <tr> <td>25</td> <td></td> </tr> </tbody> </table>		Non-Linear Limit		Value	%	25		<table border="1"> <thead> <tr> <th colspan="2">Prozone Limit</th> </tr> <tr> <th>Value</th> <th>%</th> </tr> </thead> <tbody> <tr> <td>2500.0</td> <td></td> </tr> </tbody> </table>		Prozone Limit		Value	%	2500.0		<table border="1"> <thead> <tr> <th colspan="2">User Range</th> </tr> <tr> <th>Lower</th> <th>Upper</th> </tr> </thead> <tbody> <tr> <td>3</td> <td>300</td> </tr> </tbody> </table>		User Range		Lower	Upper	3	300	<table border="1"> <thead> <tr> <th colspan="2">ValidRange</th> </tr> <tr> <th>Lower</th> <th>Upper</th> </tr> </thead> <tbody> <tr> <td>3</td> <td>300</td> </tr> </tbody> </table>		ValidRange		Lower	Upper	3	300	<table border="1"> <thead> <tr> <th colspan="2">Serum Index Limits</th> </tr> <tr> <th>Index</th> <th>Limit</th> </tr> </thead> <tbody> <tr> <td>Hemolysis</td> <td>0.0</td> </tr> <tr> <td>Icterus</td> <td>0.0</td> </tr> <tr> <td>Lipemia</td> <td>0.0</td> </tr> </tbody> </table>		Serum Index Limits		Index	Limit	Hemolysis	0.0	Icterus	0.0	Lipemia	0.0	<table border="1"> <thead> <tr> <th colspan="2">Qualitative</th> </tr> <tr> <th>On</th> <th>Off</th> </tr> </thead> <tbody> <tr> <td><input type="radio"/></td> <td><input checked="" type="radio"/></td> </tr> </tbody> </table>			Qualitative		On	Off	<input type="radio"/>	<input checked="" type="radio"/>	<table border="1"> <thead> <tr> <th colspan="2">Limit</th> <th colspan="2">Equation</th> </tr> </thead> <tbody> <tr> <td>0.0</td> <td></td> <td>none</td> <td></td> </tr> <tr> <td colspan="2"> <table border="1"> <thead> <tr> <th colspan="2">Judge Point</th> </tr> <tr> <th>Value</th> <th>%</th> </tr> </thead> <tbody> <tr> <td>0</td> <td></td> </tr> </tbody> </table> </td> <td colspan="2"></td> </tr> </tbody> </table>			Limit		Equation		0.0		none		<table border="1"> <thead> <tr> <th colspan="2">Judge Point</th> </tr> <tr> <th>Value</th> <th>%</th> </tr> </thead> <tbody> <tr> <td>0</td> <td></td> </tr> </tbody> </table>		Judge Point		Value	%	0			
Normal Range																																																																																																						
	Lower	Upper																																																																																																				
Male	0	60																																																																																																				
Female	0	60																																																																																																				
Other	0	60																																																																																																				
Reaction Slope		Absorbance Limit																																																																																																				
<input type="radio"/> Negative	<input checked="" type="radio"/> Positive	<input checked="" type="radio"/> Above	<input type="radio"/> Below																																																																																																			
<table border="1"> <thead> <tr> <th colspan="2">Non-Linear Limit</th> </tr> <tr> <th>Value</th> <th>%</th> </tr> </thead> <tbody> <tr> <td>25</td> <td></td> </tr> </tbody> </table>		Non-Linear Limit		Value	%	25		<table border="1"> <thead> <tr> <th colspan="2">Prozone Limit</th> </tr> <tr> <th>Value</th> <th>%</th> </tr> </thead> <tbody> <tr> <td>2500.0</td> <td></td> </tr> </tbody> </table>		Prozone Limit		Value	%	2500.0																																																																																								
Non-Linear Limit																																																																																																						
Value	%																																																																																																					
25																																																																																																						
Prozone Limit																																																																																																						
Value	%																																																																																																					
2500.0																																																																																																						
<table border="1"> <thead> <tr> <th colspan="2">User Range</th> </tr> <tr> <th>Lower</th> <th>Upper</th> </tr> </thead> <tbody> <tr> <td>3</td> <td>300</td> </tr> </tbody> </table>		User Range		Lower	Upper	3	300	<table border="1"> <thead> <tr> <th colspan="2">ValidRange</th> </tr> <tr> <th>Lower</th> <th>Upper</th> </tr> </thead> <tbody> <tr> <td>3</td> <td>300</td> </tr> </tbody> </table>		ValidRange		Lower	Upper	3	300	<table border="1"> <thead> <tr> <th colspan="2">Serum Index Limits</th> </tr> <tr> <th>Index</th> <th>Limit</th> </tr> </thead> <tbody> <tr> <td>Hemolysis</td> <td>0.0</td> </tr> <tr> <td>Icterus</td> <td>0.0</td> </tr> <tr> <td>Lipemia</td> <td>0.0</td> </tr> </tbody> </table>		Serum Index Limits		Index	Limit	Hemolysis	0.0	Icterus	0.0	Lipemia	0.0																																																																											
User Range																																																																																																						
Lower	Upper																																																																																																					
3	300																																																																																																					
ValidRange																																																																																																						
Lower	Upper																																																																																																					
3	300																																																																																																					
Serum Index Limits																																																																																																						
Index	Limit																																																																																																					
Hemolysis	0.0																																																																																																					
Icterus	0.0																																																																																																					
Lipemia	0.0																																																																																																					
<table border="1"> <thead> <tr> <th colspan="2">Qualitative</th> </tr> <tr> <th>On</th> <th>Off</th> </tr> </thead> <tbody> <tr> <td><input type="radio"/></td> <td><input checked="" type="radio"/></td> </tr> </tbody> </table>			Qualitative		On	Off	<input type="radio"/>	<input checked="" type="radio"/>	<table border="1"> <thead> <tr> <th colspan="2">Limit</th> <th colspan="2">Equation</th> </tr> </thead> <tbody> <tr> <td>0.0</td> <td></td> <td>none</td> <td></td> </tr> <tr> <td colspan="2"> <table border="1"> <thead> <tr> <th colspan="2">Judge Point</th> </tr> <tr> <th>Value</th> <th>%</th> </tr> </thead> <tbody> <tr> <td>0</td> <td></td> </tr> </tbody> </table> </td> <td colspan="2"></td> </tr> </tbody> </table>			Limit		Equation		0.0		none		<table border="1"> <thead> <tr> <th colspan="2">Judge Point</th> </tr> <tr> <th>Value</th> <th>%</th> </tr> </thead> <tbody> <tr> <td>0</td> <td></td> </tr> </tbody> </table>		Judge Point		Value	%	0																																																																												
Qualitative																																																																																																						
On	Off																																																																																																					
<input type="radio"/>	<input checked="" type="radio"/>																																																																																																					
Limit		Equation																																																																																																				
0.0		none																																																																																																				
<table border="1"> <thead> <tr> <th colspan="2">Judge Point</th> </tr> <tr> <th>Value</th> <th>%</th> </tr> </thead> <tbody> <tr> <td>0</td> <td></td> </tr> </tbody> </table>		Judge Point		Value	%	0																																																																																																
Judge Point																																																																																																						
Value	%																																																																																																					
0																																																																																																						

Test Name		Test Code		Measure		Sample Reagent		Ranges Limits		Calibration																																																																																											
NN	Lipase	IL	Lps -d																																																																																																		
Calibration				Condition				Calibration/R-Blank Limit																																																																																													
<table border="1"> <thead> <tr> <th colspan="2">Method</th> <th colspan="2">Curve Type</th> </tr> </thead> <tbody> <tr> <td>1-point</td> <td></td> <td>Linear</td> <td></td> </tr> <tr> <td colspan="2"> <table border="1"> <thead> <tr> <th colspan="2">Repeats</th> </tr> <tr> <th>Value</th> <th>%</th> </tr> </thead> <tbody> <tr> <td>2</td> <td></td> </tr> </tbody> </table> </td> <td colspan="2"> <table border="1"> <thead> <tr> <th colspan="2">Stability</th> </tr> <tr> <th>Value</th> <th>days</th> </tr> </thead> <tbody> <tr> <td></td> <td></td> </tr> </tbody> </table> </td> </tr> <tr> <td colspan="4"> <table border="1"> <thead> <tr> <th colspan="2">Calibrator</th> </tr> <tr> <th>Calibrator</th> <th>Conc</th> </tr> </thead> <tbody> <tr> <td>1</td> <td>TruCalU **</td> </tr> <tr> <td>2</td> <td></td> </tr> <tr> <td>3</td> <td></td> </tr> <tr> <td>4</td> <td></td> </tr> <tr> <td>5</td> <td></td> </tr> </tbody> </table> </td> <td colspan="4"> <table border="1"> <thead> <tr> <th colspan="2">Condition</th> </tr> </thead> <tbody> <tr> <td><input checked="" type="checkbox"/></td> <td>Reagent Blank</td> </tr> <tr> <td><input type="checkbox"/></td> <td>Auto Reagent Blank by Bottle</td> </tr> <tr> <td><input type="checkbox"/></td> <td>Ask for calibration when reagent lot changes</td> </tr> </tbody> </table> </td> <td colspan="4"> <table border="1"> <thead> <tr> <th colspan="2">R-Blank Limit</th> </tr> <tr> <th>Value</th> <th>mAbs</th> </tr> </thead> <tbody> <tr> <td>2500.0</td> <td></td> </tr> <tr> <td colspan="2"> <table border="1"> <thead> <tr> <th colspan="2">Cal Reps Range</th> </tr> <tr> <th>Value</th> <th>%</th> </tr> </thead> <tbody> <tr> <td>15.0</td> <td></td> </tr> </tbody> </table> </td> </tr> <tr> <td colspan="2"> <table border="1"> <thead> <tr> <th colspan="2">Min Cal Reps</th> </tr> <tr> <th>Value</th> <th>mAbs</th> </tr> </thead> <tbody> <tr> <td>50.0</td> <td></td> </tr> </tbody> </table> </td> <td colspan="2"> <table border="1"> <thead> <tr> <th colspan="2">Factor Change</th> </tr> <tr> <th>Value</th> <th>%</th> </tr> </thead> <tbody> <tr> <td>20</td> <td></td> </tr> </tbody> </table> </td> </tr> <tr> <td colspan="2"> <table border="1"> <thead> <tr> <th colspan="2">M-Point Curve Fit</th> </tr> <tr> <th>Value</th> <th>%</th> </tr> </thead> <tbody> <tr> <td>0.0</td> <td></td> </tr> </tbody> </table> </td> <td colspan="2"></td> </tr> </tbody> </table> </td> </tr> </tbody> </table>				Method		Curve Type		1-point		Linear		<table border="1"> <thead> <tr> <th colspan="2">Repeats</th> </tr> <tr> <th>Value</th> <th>%</th> </tr> </thead> <tbody> <tr> <td>2</td> <td></td> </tr> </tbody> </table>		Repeats		Value	%	2		<table border="1"> <thead> <tr> <th colspan="2">Stability</th> </tr> <tr> <th>Value</th> <th>days</th> </tr> </thead> <tbody> <tr> <td></td> <td></td> </tr> </tbody> </table>		Stability		Value	days			<table border="1"> <thead> <tr> <th colspan="2">Calibrator</th> </tr> <tr> <th>Calibrator</th> <th>Conc</th> </tr> </thead> <tbody> <tr> <td>1</td> <td>TruCalU **</td> </tr> <tr> <td>2</td> <td></td> </tr> <tr> <td>3</td> <td></td> </tr> <tr> <td>4</td> <td></td> </tr> <tr> <td>5</td> <td></td> </tr> </tbody> </table>				Calibrator		Calibrator	Conc	1	TruCalU **	2		3		4		5		<table border="1"> <thead> <tr> <th colspan="2">Condition</th> </tr> </thead> <tbody> <tr> <td><input checked="" type="checkbox"/></td> <td>Reagent Blank</td> </tr> <tr> <td><input type="checkbox"/></td> <td>Auto Reagent Blank by Bottle</td> </tr> <tr> <td><input type="checkbox"/></td> <td>Ask for calibration when reagent lot changes</td> </tr> </tbody> </table>				Condition		<input checked="" type="checkbox"/>	Reagent Blank	<input type="checkbox"/>	Auto Reagent Blank by Bottle	<input type="checkbox"/>	Ask for calibration when reagent lot changes	<table border="1"> <thead> <tr> <th colspan="2">R-Blank Limit</th> </tr> <tr> <th>Value</th> <th>mAbs</th> </tr> </thead> <tbody> <tr> <td>2500.0</td> <td></td> </tr> <tr> <td colspan="2"> <table border="1"> <thead> <tr> <th colspan="2">Cal Reps Range</th> </tr> <tr> <th>Value</th> <th>%</th> </tr> </thead> <tbody> <tr> <td>15.0</td> <td></td> </tr> </tbody> </table> </td> </tr> <tr> <td colspan="2"> <table border="1"> <thead> <tr> <th colspan="2">Min Cal Reps</th> </tr> <tr> <th>Value</th> <th>mAbs</th> </tr> </thead> <tbody> <tr> <td>50.0</td> <td></td> </tr> </tbody> </table> </td> <td colspan="2"> <table border="1"> <thead> <tr> <th colspan="2">Factor Change</th> </tr> <tr> <th>Value</th> <th>%</th> </tr> </thead> <tbody> <tr> <td>20</td> <td></td> </tr> </tbody> </table> </td> </tr> <tr> <td colspan="2"> <table border="1"> <thead> <tr> <th colspan="2">M-Point Curve Fit</th> </tr> <tr> <th>Value</th> <th>%</th> </tr> </thead> <tbody> <tr> <td>0.0</td> <td></td> </tr> </tbody> </table> </td> <td colspan="2"></td> </tr> </tbody> </table>				R-Blank Limit		Value	mAbs	2500.0		<table border="1"> <thead> <tr> <th colspan="2">Cal Reps Range</th> </tr> <tr> <th>Value</th> <th>%</th> </tr> </thead> <tbody> <tr> <td>15.0</td> <td></td> </tr> </tbody> </table>		Cal Reps Range		Value	%	15.0		<table border="1"> <thead> <tr> <th colspan="2">Min Cal Reps</th> </tr> <tr> <th>Value</th> <th>mAbs</th> </tr> </thead> <tbody> <tr> <td>50.0</td> <td></td> </tr> </tbody> </table>		Min Cal Reps		Value	mAbs	50.0		<table border="1"> <thead> <tr> <th colspan="2">Factor Change</th> </tr> <tr> <th>Value</th> <th>%</th> </tr> </thead> <tbody> <tr> <td>20</td> <td></td> </tr> </tbody> </table>		Factor Change		Value	%	20		<table border="1"> <thead> <tr> <th colspan="2">M-Point Curve Fit</th> </tr> <tr> <th>Value</th> <th>%</th> </tr> </thead> <tbody> <tr> <td>0.0</td> <td></td> </tr> </tbody> </table>		M-Point Curve Fit		Value	%	0.0			
Method		Curve Type																																																																																																			
1-point		Linear																																																																																																			
<table border="1"> <thead> <tr> <th colspan="2">Repeats</th> </tr> <tr> <th>Value</th> <th>%</th> </tr> </thead> <tbody> <tr> <td>2</td> <td></td> </tr> </tbody> </table>		Repeats		Value	%	2		<table border="1"> <thead> <tr> <th colspan="2">Stability</th> </tr> <tr> <th>Value</th> <th>days</th> </tr> </thead> <tbody> <tr> <td></td> <td></td> </tr> </tbody> </table>		Stability		Value	days																																																																																								
Repeats																																																																																																					
Value	%																																																																																																				
2																																																																																																					
Stability																																																																																																					
Value	days																																																																																																				
<table border="1"> <thead> <tr> <th colspan="2">Calibrator</th> </tr> <tr> <th>Calibrator</th> <th>Conc</th> </tr> </thead> <tbody> <tr> <td>1</td> <td>TruCalU **</td> </tr> <tr> <td>2</td> <td></td> </tr> <tr> <td>3</td> <td></td> </tr> <tr> <td>4</td> <td></td> </tr> <tr> <td>5</td> <td></td> </tr> </tbody> </table>				Calibrator		Calibrator	Conc	1	TruCalU **	2		3		4		5		<table border="1"> <thead> <tr> <th colspan="2">Condition</th> </tr> </thead> <tbody> <tr> <td><input checked="" type="checkbox"/></td> <td>Reagent Blank</td> </tr> <tr> <td><input type="checkbox"/></td> <td>Auto Reagent Blank by Bottle</td> </tr> <tr> <td><input type="checkbox"/></td> <td>Ask for calibration when reagent lot changes</td> </tr> </tbody> </table>				Condition		<input checked="" type="checkbox"/>	Reagent Blank	<input type="checkbox"/>	Auto Reagent Blank by Bottle	<input type="checkbox"/>	Ask for calibration when reagent lot changes	<table border="1"> <thead> <tr> <th colspan="2">R-Blank Limit</th> </tr> <tr> <th>Value</th> <th>mAbs</th> </tr> </thead> <tbody> <tr> <td>2500.0</td> <td></td> </tr> <tr> <td colspan="2"> <table border="1"> <thead> <tr> <th colspan="2">Cal Reps Range</th> </tr> <tr> <th>Value</th> <th>%</th> </tr> </thead> <tbody> <tr> <td>15.0</td> <td></td> </tr> </tbody> </table> </td> </tr> <tr> <td colspan="2"> <table border="1"> <thead> <tr> <th colspan="2">Min Cal Reps</th> </tr> <tr> <th>Value</th> <th>mAbs</th> </tr> </thead> <tbody> <tr> <td>50.0</td> <td></td> </tr> </tbody> </table> </td> <td colspan="2"> <table border="1"> <thead> <tr> <th colspan="2">Factor Change</th> </tr> <tr> <th>Value</th> <th>%</th> </tr> </thead> <tbody> <tr> <td>20</td> <td></td> </tr> </tbody> </table> </td> </tr> <tr> <td colspan="2"> <table border="1"> <thead> <tr> <th colspan="2">M-Point Curve Fit</th> </tr> <tr> <th>Value</th> <th>%</th> </tr> </thead> <tbody> <tr> <td>0.0</td> <td></td> </tr> </tbody> </table> </td> <td colspan="2"></td> </tr> </tbody> </table>				R-Blank Limit		Value	mAbs	2500.0		<table border="1"> <thead> <tr> <th colspan="2">Cal Reps Range</th> </tr> <tr> <th>Value</th> <th>%</th> </tr> </thead> <tbody> <tr> <td>15.0</td> <td></td> </tr> </tbody> </table>		Cal Reps Range		Value	%	15.0		<table border="1"> <thead> <tr> <th colspan="2">Min Cal Reps</th> </tr> <tr> <th>Value</th> <th>mAbs</th> </tr> </thead> <tbody> <tr> <td>50.0</td> <td></td> </tr> </tbody> </table>		Min Cal Reps		Value	mAbs	50.0		<table border="1"> <thead> <tr> <th colspan="2">Factor Change</th> </tr> <tr> <th>Value</th> <th>%</th> </tr> </thead> <tbody> <tr> <td>20</td> <td></td> </tr> </tbody> </table>		Factor Change		Value	%	20		<table border="1"> <thead> <tr> <th colspan="2">M-Point Curve Fit</th> </tr> <tr> <th>Value</th> <th>%</th> </tr> </thead> <tbody> <tr> <td>0.0</td> <td></td> </tr> </tbody> </table>		M-Point Curve Fit		Value	%	0.0																															
Calibrator																																																																																																					
Calibrator	Conc																																																																																																				
1	TruCalU **																																																																																																				
2																																																																																																					
3																																																																																																					
4																																																																																																					
5																																																																																																					
Condition																																																																																																					
<input checked="" type="checkbox"/>	Reagent Blank																																																																																																				
<input type="checkbox"/>	Auto Reagent Blank by Bottle																																																																																																				
<input type="checkbox"/>	Ask for calibration when reagent lot changes																																																																																																				
R-Blank Limit																																																																																																					
Value	mAbs																																																																																																				
2500.0																																																																																																					
<table border="1"> <thead> <tr> <th colspan="2">Cal Reps Range</th> </tr> <tr> <th>Value</th> <th>%</th> </tr> </thead> <tbody> <tr> <td>15.0</td> <td></td> </tr> </tbody> </table>		Cal Reps Range		Value	%	15.0																																																																																															
Cal Reps Range																																																																																																					
Value	%																																																																																																				
15.0																																																																																																					
<table border="1"> <thead> <tr> <th colspan="2">Min Cal Reps</th> </tr> <tr> <th>Value</th> <th>mAbs</th> </tr> </thead> <tbody> <tr> <td>50.0</td> <td></td> </tr> </tbody> </table>		Min Cal Reps		Value	mAbs	50.0		<table border="1"> <thead> <tr> <th colspan="2">Factor Change</th> </tr> <tr> <th>Value</th> <th>%</th> </tr> </thead> <tbody> <tr> <td>20</td> <td></td> </tr> </tbody> </table>		Factor Change		Value	%	20																																																																																							
Min Cal Reps																																																																																																					
Value	mAbs																																																																																																				
50.0																																																																																																					
Factor Change																																																																																																					
Value	%																																																																																																				
20																																																																																																					
<table border="1"> <thead> <tr> <th colspan="2">M-Point Curve Fit</th> </tr> <tr> <th>Value</th> <th>%</th> </tr> </thead> <tbody> <tr> <td>0.0</td> <td></td> </tr> </tbody> </table>		M-Point Curve Fit		Value	%	0.0																																																																																															
M-Point Curve Fit																																																																																																					
Value	%																																																																																																				
0.0																																																																																																					

* -Рекомендуемый параметр.

** - Вводится из паспорта калибратору.

Тип калибровки: Линейная (калибратор TruCalU) или по фактору (см инструкцию).
Контроль по TruLab N и TruLab P.