

## Параметры для ввода в программу анализатора iLab-Taurus

Test Name: NN CI Test Code: IL\_ CI-d Measure Sample Reagent Ranges Limits Calibration

Sample Type: Serum

Reaction Cycle: Standard Extended

Reporting Unit: ммоль/л \* Decimal Points: 1

User Define: Conv.Factor: 0.000

Methodology: Type: End Point Rate Measuring Point: 10 Photometric: 1 Wavelength 2 Wavelength Primary: 450 Secondary: 600

Correction Constant: Slope: 1.000 Intercept: 0.000

Test Name: NN CI Test Code: IL\_ CI-d Measure Sample Reagent Ranges Limits Calibration

Sample Volume: Sampling Condition:

Sample	Dilution	Diluent	Condition No.	1	2	3	4
1	2.0	0.0	First Run				
2	2.0	20.0	Samp.Vol.Reduction				
3	4.0	0.0	Below N-Range				
4			Above N-Range				
			Panic L				
			Panic H				
			User Range L				
			User Range H				
			Noise				
			Prozone				
			HIGH!				
			ABS!				

Diluent: Diluent Warning Limit: 0 Tests

Reagent Volume: R1: CI-d R2:

	Volume	Diluent Vol.	stirring	Warning Limit(tests)	Stability(days)
R1	200	0.0		20	No Control
R2					

# Хлориды

Thiocyanate (тиоцианатный метод)

Ranges	Limits																																																																					
<table><thead><tr><th colspan="3">Normal Range</th></tr><tr><th></th><th>Lower</th><th>Upper</th></tr></thead><tbody><tr><td>Male</td><td>95</td><td>105</td></tr><tr><td>Female</td><td>95</td><td>105</td></tr><tr><td>Other</td><td>95</td><td>105</td></tr></tbody></table> <table><thead><tr><th colspan="2">User Range</th></tr><tr><th>Lower</th><th>Upper</th></tr></thead><tbody><tr><td>1</td><td>150</td></tr></tbody></table> <table><thead><tr><th colspan="2">ValidRange</th></tr><tr><th>Lower</th><th>Upper</th></tr></thead><tbody><tr><td>1</td><td>150</td></tr></tbody></table>	Normal Range				Lower	Upper	Male	95	105	Female	95	105	Other	95	105	User Range		Lower	Upper	1	150	ValidRange		Lower	Upper	1	150	<table><thead><tr><th>Reaction Slope</th><th>Absorbance Limit</th></tr></thead><tbody><tr><td><input type="radio"/> Negative <input checked="" type="radio"/> Positive</td><td><input checked="" type="radio"/> Above <input type="radio"/> Below</td></tr><tr><td><table><thead><tr><th colspan="2">Non-Linear Limit</th></tr><tr><td>0</td><td>%</td></tr></thead></table></td><td><table><tr><td>2500.0</td><td>mAbs</td></tr></table></td></tr><tr><td colspan="2"><table><thead><tr><th colspan="2">Serum Index Limits</th></tr><tr><th></th><th>Limit</th></tr></thead><tbody><tr><td>Hemolysis</td><td>0.0</td></tr><tr><td>Icterus</td><td>0.0</td></tr><tr><td>Lipemia</td><td>0.0</td></tr></tbody></table></td></tr><tr><td colspan="2"><table><thead><tr><th colspan="2">Prozone Limit</th></tr><tr><td><input type="radio"/> Above <input checked="" type="radio"/> Below</td><td></td></tr><tr><td><table><thead><tr><th>Limit</th><th>Equation</th></tr></thead><tbody><tr><td>0.0</td><td>none</td></tr></tbody></table></td><td></td></tr><tr><td><table><thead><tr><th colspan="2">Judge Point</th></tr><tr><td>0</td><td></td></tr></thead></table></td><td></td></tr></thead></table></td></tr></tbody></table>	Reaction Slope	Absorbance Limit	<input type="radio"/> Negative <input checked="" type="radio"/> Positive	<input checked="" type="radio"/> Above <input type="radio"/> Below	<table><thead><tr><th colspan="2">Non-Linear Limit</th></tr><tr><td>0</td><td>%</td></tr></thead></table>	Non-Linear Limit		0	%	<table><tr><td>2500.0</td><td>mAbs</td></tr></table>	2500.0	mAbs	<table><thead><tr><th colspan="2">Serum Index Limits</th></tr><tr><th></th><th>Limit</th></tr></thead><tbody><tr><td>Hemolysis</td><td>0.0</td></tr><tr><td>Icterus</td><td>0.0</td></tr><tr><td>Lipemia</td><td>0.0</td></tr></tbody></table>		Serum Index Limits			Limit	Hemolysis	0.0	Icterus	0.0	Lipemia	0.0	<table><thead><tr><th colspan="2">Prozone Limit</th></tr><tr><td><input type="radio"/> Above <input checked="" type="radio"/> Below</td><td></td></tr><tr><td><table><thead><tr><th>Limit</th><th>Equation</th></tr></thead><tbody><tr><td>0.0</td><td>none</td></tr></tbody></table></td><td></td></tr><tr><td><table><thead><tr><th colspan="2">Judge Point</th></tr><tr><td>0</td><td></td></tr></thead></table></td><td></td></tr></thead></table>		Prozone Limit		<input type="radio"/> Above <input checked="" type="radio"/> Below		<table><thead><tr><th>Limit</th><th>Equation</th></tr></thead><tbody><tr><td>0.0</td><td>none</td></tr></tbody></table>	Limit	Equation	0.0	none		<table><thead><tr><th colspan="2">Judge Point</th></tr><tr><td>0</td><td></td></tr></thead></table>	Judge Point		0		
Normal Range																																																																						
	Lower	Upper																																																																				
Male	95	105																																																																				
Female	95	105																																																																				
Other	95	105																																																																				
User Range																																																																						
Lower	Upper																																																																					
1	150																																																																					
ValidRange																																																																						
Lower	Upper																																																																					
1	150																																																																					
Reaction Slope	Absorbance Limit																																																																					
<input type="radio"/> Negative <input checked="" type="radio"/> Positive	<input checked="" type="radio"/> Above <input type="radio"/> Below																																																																					
<table><thead><tr><th colspan="2">Non-Linear Limit</th></tr><tr><td>0</td><td>%</td></tr></thead></table>	Non-Linear Limit		0	%	<table><tr><td>2500.0</td><td>mAbs</td></tr></table>	2500.0	mAbs																																																															
Non-Linear Limit																																																																						
0	%																																																																					
2500.0	mAbs																																																																					
<table><thead><tr><th colspan="2">Serum Index Limits</th></tr><tr><th></th><th>Limit</th></tr></thead><tbody><tr><td>Hemolysis</td><td>0.0</td></tr><tr><td>Icterus</td><td>0.0</td></tr><tr><td>Lipemia</td><td>0.0</td></tr></tbody></table>		Serum Index Limits			Limit	Hemolysis	0.0	Icterus	0.0	Lipemia	0.0																																																											
Serum Index Limits																																																																						
	Limit																																																																					
Hemolysis	0.0																																																																					
Icterus	0.0																																																																					
Lipemia	0.0																																																																					
<table><thead><tr><th colspan="2">Prozone Limit</th></tr><tr><td><input type="radio"/> Above <input checked="" type="radio"/> Below</td><td></td></tr><tr><td><table><thead><tr><th>Limit</th><th>Equation</th></tr></thead><tbody><tr><td>0.0</td><td>none</td></tr></tbody></table></td><td></td></tr><tr><td><table><thead><tr><th colspan="2">Judge Point</th></tr><tr><td>0</td><td></td></tr></thead></table></td><td></td></tr></thead></table>		Prozone Limit		<input type="radio"/> Above <input checked="" type="radio"/> Below		<table><thead><tr><th>Limit</th><th>Equation</th></tr></thead><tbody><tr><td>0.0</td><td>none</td></tr></tbody></table>	Limit	Equation	0.0	none		<table><thead><tr><th colspan="2">Judge Point</th></tr><tr><td>0</td><td></td></tr></thead></table>	Judge Point		0																																																							
Prozone Limit																																																																						
<input type="radio"/> Above <input checked="" type="radio"/> Below																																																																						
<table><thead><tr><th>Limit</th><th>Equation</th></tr></thead><tbody><tr><td>0.0</td><td>none</td></tr></tbody></table>	Limit	Equation	0.0	none																																																																		
Limit	Equation																																																																					
0.0	none																																																																					
<table><thead><tr><th colspan="2">Judge Point</th></tr><tr><td>0</td><td></td></tr></thead></table>	Judge Point		0																																																																			
Judge Point																																																																						
0																																																																						
<table><thead><tr><th colspan="2">Qualitative</th></tr><tr><td><input type="radio"/> On <input checked="" type="radio"/> Off</td><td></td></tr></thead></table>	Qualitative		<input type="radio"/> On <input checked="" type="radio"/> Off																																																																			
Qualitative																																																																						
<input type="radio"/> On <input checked="" type="radio"/> Off																																																																						

Test Name	Test Code	Measure	Sample Reagent	Ranges Limits	Calibration																																																				
NN	CI	IL_	CI-d																																																						
Calibration		Condition		Calibration/R-Blank Limit																																																					
<table><thead><tr><th>Method</th><th>Curve Type</th></tr></thead><tbody><tr><td>1-point</td><td>Linear</td></tr></tbody></table> <table><thead><tr><th>Repeats</th><th>Stability</th></tr></thead><tbody><tr><td>2</td><td>days</td></tr></tbody></table> <table><thead><tr><th>No</th></tr></thead><tbody><tr><td></td></tr></tbody></table> <table><thead><tr><th colspan="2">Calibrator</th></tr><tr><th>Calibrator</th><th>Conc</th></tr></thead><tbody><tr><td>1</td><td>TruCalU **</td></tr><tr><td>2</td><td></td></tr><tr><td>3</td><td></td></tr><tr><td>4</td><td></td></tr><tr><td>5</td><td></td></tr></tbody></table>		Method	Curve Type	1-point	Linear	Repeats	Stability	2	days	No		Calibrator		Calibrator	Conc	1	TruCalU **	2		3		4		5		<input checked="" type="checkbox"/> Reagent Blank <input type="checkbox"/> Auto Reagent Blank by Bottle <input type="checkbox"/> Ask for calibration when reagent lot changes		<table><thead><tr><th colspan="2">R-Blank Limit</th></tr><tr><td>2500.0</td><td>mAbs</td></tr></thead><tbody><tr><td colspan="2"><table><thead><tr><th colspan="2">Cal Reps Range</th></tr><tr><td>20.0</td><td>%</td></tr></thead></table></td></tr><tr><td colspan="2"><table><thead><tr><th colspan="2">Min Cal Reps</th></tr><tr><td>50.0</td><td>mAbs</td></tr></thead></table></td></tr><tr><td colspan="2"><table><thead><tr><th colspan="2">Factor Change</th></tr><tr><td>20</td><td>%</td></tr></thead></table></td></tr><tr><td colspan="2"><table><thead><tr><th colspan="2">M-Point Curve Fit</th></tr><tr><td>0.0</td><td>%</td></tr></thead></table></td></tr></tbody></table>		R-Blank Limit		2500.0	mAbs	<table><thead><tr><th colspan="2">Cal Reps Range</th></tr><tr><td>20.0</td><td>%</td></tr></thead></table>		Cal Reps Range		20.0	%	<table><thead><tr><th colspan="2">Min Cal Reps</th></tr><tr><td>50.0</td><td>mAbs</td></tr></thead></table>		Min Cal Reps		50.0	mAbs	<table><thead><tr><th colspan="2">Factor Change</th></tr><tr><td>20</td><td>%</td></tr></thead></table>		Factor Change		20	%	<table><thead><tr><th colspan="2">M-Point Curve Fit</th></tr><tr><td>0.0</td><td>%</td></tr></thead></table>		M-Point Curve Fit		0.0	%
Method	Curve Type																																																								
1-point	Linear																																																								
Repeats	Stability																																																								
2	days																																																								
No																																																									
Calibrator																																																									
Calibrator	Conc																																																								
1	TruCalU **																																																								
2																																																									
3																																																									
4																																																									
5																																																									
R-Blank Limit																																																									
2500.0	mAbs																																																								
<table><thead><tr><th colspan="2">Cal Reps Range</th></tr><tr><td>20.0</td><td>%</td></tr></thead></table>		Cal Reps Range		20.0	%																																																				
Cal Reps Range																																																									
20.0	%																																																								
<table><thead><tr><th colspan="2">Min Cal Reps</th></tr><tr><td>50.0</td><td>mAbs</td></tr></thead></table>		Min Cal Reps		50.0	mAbs																																																				
Min Cal Reps																																																									
50.0	mAbs																																																								
<table><thead><tr><th colspan="2">Factor Change</th></tr><tr><td>20</td><td>%</td></tr></thead></table>		Factor Change		20	%																																																				
Factor Change																																																									
20	%																																																								
<table><thead><tr><th colspan="2">M-Point Curve Fit</th></tr><tr><td>0.0</td><td>%</td></tr></thead></table>		M-Point Curve Fit		0.0	%																																																				
M-Point Curve Fit																																																									
0.0	%																																																								

# Хлориды

*Thiocyanate (тиоцианатный метод)*



\* -Рекомендуемый параметр.

\*\* - Вводится из паспорта калибратору.

Тип калибровки: Линейная (калибратор TruCalU) или по стандарту из набора.

Контроль по TruLab N и TruLab P.