

## Параметры для ввода в программу анализатора iLab-Taurus

Test Name: NN CK Test Code: IL\_ CK-d

Measure Sample Reagent Ranges Limits Calibration

Sample Type: Serum

Reaction Cycle: Standard Extended

Reporting Unit: Ед/л \* Decimal Points: 1

User Define: Conv.Factor: 0.000

Methodology: Type: End Point Rate

Measuring Point: 21 - 33

Photometric: 1 Wavelength 2 Wavelength

Primary: 340 Secondary: 546

Correction Constant: Slope: 1.000 Intercept: 0.000

Test Name: NN CK Test Code: IL\_ CK-d

Measure Sample Reagent Ranges Limits Calibration

Sample Volume: Sampling Condition:

Sample	Dilution	Diluent	Condition No.	1	2	3	4
1	7.0	0.0	First Run				
2	3.0	0.0	Samp.Vol.Reduction				
3	14.0	0.0	Below N-Range				
4	2.0	20	Above N-Range				
			Panic L				
			Panic H				
			User Range L				
			User Range H				
			Noise				
			Prozone				
			HIGH!				
			ABS!				

Diluent: Diluent Warning Limit: 0 Tests

Reagent Volume: R1: CK-d1 R2: CK-d2

	Volume	Diluent Vol.	stirring	Warning Limit(tests)	Stability(days)
R1	140	0.0		20	No Control
R2	35	0.0		20	No Control

Ranges		Limits																																																																																												
<table border="1"> <thead> <tr> <th colspan="3">Normal Range</th> </tr> <tr> <th></th> <th>Lower</th> <th>Upper</th> </tr> </thead> <tbody> <tr> <td>Male</td> <td>0</td> <td>171</td> </tr> <tr> <td>Female</td> <td>0</td> <td>145</td> </tr> <tr> <td>Other</td> <td>0</td> <td>171</td> </tr> </tbody> </table>		Normal Range				Lower	Upper	Male	0	171	Female	0	145	Other	0	171	<table border="1"> <thead> <tr> <th colspan="2">Reaction Slope</th> <th colspan="2">Absorbance Limit</th> </tr> </thead> <tbody> <tr> <td colspan="2"> <input type="radio"/> Negative <input checked="" type="radio"/> Positive                 </td> <td colspan="2"> <input checked="" type="radio"/> Above <input type="radio"/> Below             </td> </tr> <tr> <td colspan="2"> <table border="1"> <thead> <tr> <th colspan="2">Non-Linear Limit</th> </tr> <tr> <th>Value</th> <th>%</th> </tr> </thead> <tbody> <tr> <td>25</td> <td></td> </tr> </tbody> </table> </td> <td colspan="2"> <table border="1"> <thead> <tr> <th colspan="2">Prozone Limit</th> </tr> <tr> <th>Value</th> <th>%</th> </tr> </thead> <tbody> <tr> <td>2500.0</td> <td></td> </tr> </tbody> </table> </td> </tr> <tr> <td colspan="2"> <table border="1"> <thead> <tr> <th colspan="2">User Range</th> </tr> <tr> <th>Lower</th> <th>Upper</th> </tr> </thead> <tbody> <tr> <td>1</td> <td>1100</td> </tr> </tbody> </table> </td> <td colspan="2"> <table border="1"> <thead> <tr> <th colspan="2">Serum Index Limits</th> </tr> <tr> <th></th> <th>Limit</th> </tr> </thead> <tbody> <tr> <td>Hemolysis</td> <td>0.0</td> </tr> <tr> <td>Icterus</td> <td>0.0</td> </tr> <tr> <td>Lipemia</td> <td>0.0</td> </tr> </tbody> </table> </td> </tr> <tr> <td colspan="2"> <table border="1"> <thead> <tr> <th colspan="2">ValidRange</th> </tr> <tr> <th>Lower</th> <th>Upper</th> </tr> </thead> <tbody> <tr> <td>1</td> <td>1100</td> </tr> </tbody> </table> </td> <td colspan="2"> <table border="1"> <thead> <tr> <th colspan="2">Equation</th> </tr> <tr> <th>Limit</th> <th>Equation</th> </tr> </thead> <tbody> <tr> <td>0.0</td> <td>none</td> </tr> </tbody> </table> </td> </tr> <tr> <td colspan="2"> <table border="1"> <thead> <tr> <th colspan="2">Qualitative</th> </tr> <tr> <th>On</th> <th>Off</th> </tr> </thead> <tbody> <tr> <td><input type="radio"/></td> <td><input checked="" type="radio"/></td> </tr> </tbody> </table> </td> <td colspan="2"> <table border="1"> <thead> <tr> <th colspan="2">Judge Point</th> </tr> <tr> <th>Value</th> <th>%</th> </tr> </thead> <tbody> <tr> <td>0</td> <td></td> </tr> </tbody> </table> </td> </tr> </tbody> </table>		Reaction Slope		Absorbance Limit		<input type="radio"/> Negative <input checked="" type="radio"/> Positive		<input checked="" type="radio"/> Above <input type="radio"/> Below		<table border="1"> <thead> <tr> <th colspan="2">Non-Linear Limit</th> </tr> <tr> <th>Value</th> <th>%</th> </tr> </thead> <tbody> <tr> <td>25</td> <td></td> </tr> </tbody> </table>		Non-Linear Limit		Value	%	25		<table border="1"> <thead> <tr> <th colspan="2">Prozone Limit</th> </tr> <tr> <th>Value</th> <th>%</th> </tr> </thead> <tbody> <tr> <td>2500.0</td> <td></td> </tr> </tbody> </table>		Prozone Limit		Value	%	2500.0		<table border="1"> <thead> <tr> <th colspan="2">User Range</th> </tr> <tr> <th>Lower</th> <th>Upper</th> </tr> </thead> <tbody> <tr> <td>1</td> <td>1100</td> </tr> </tbody> </table>		User Range		Lower	Upper	1	1100	<table border="1"> <thead> <tr> <th colspan="2">Serum Index Limits</th> </tr> <tr> <th></th> <th>Limit</th> </tr> </thead> <tbody> <tr> <td>Hemolysis</td> <td>0.0</td> </tr> <tr> <td>Icterus</td> <td>0.0</td> </tr> <tr> <td>Lipemia</td> <td>0.0</td> </tr> </tbody> </table>		Serum Index Limits			Limit	Hemolysis	0.0	Icterus	0.0	Lipemia	0.0	<table border="1"> <thead> <tr> <th colspan="2">ValidRange</th> </tr> <tr> <th>Lower</th> <th>Upper</th> </tr> </thead> <tbody> <tr> <td>1</td> <td>1100</td> </tr> </tbody> </table>		ValidRange		Lower	Upper	1	1100	<table border="1"> <thead> <tr> <th colspan="2">Equation</th> </tr> <tr> <th>Limit</th> <th>Equation</th> </tr> </thead> <tbody> <tr> <td>0.0</td> <td>none</td> </tr> </tbody> </table>		Equation		Limit	Equation	0.0	none	<table border="1"> <thead> <tr> <th colspan="2">Qualitative</th> </tr> <tr> <th>On</th> <th>Off</th> </tr> </thead> <tbody> <tr> <td><input type="radio"/></td> <td><input checked="" type="radio"/></td> </tr> </tbody> </table>		Qualitative		On	Off	<input type="radio"/>	<input checked="" type="radio"/>	<table border="1"> <thead> <tr> <th colspan="2">Judge Point</th> </tr> <tr> <th>Value</th> <th>%</th> </tr> </thead> <tbody> <tr> <td>0</td> <td></td> </tr> </tbody> </table>		Judge Point		Value	%	0	
Normal Range																																																																																														
	Lower	Upper																																																																																												
Male	0	171																																																																																												
Female	0	145																																																																																												
Other	0	171																																																																																												
Reaction Slope		Absorbance Limit																																																																																												
<input type="radio"/> Negative <input checked="" type="radio"/> Positive		<input checked="" type="radio"/> Above <input type="radio"/> Below																																																																																												
<table border="1"> <thead> <tr> <th colspan="2">Non-Linear Limit</th> </tr> <tr> <th>Value</th> <th>%</th> </tr> </thead> <tbody> <tr> <td>25</td> <td></td> </tr> </tbody> </table>		Non-Linear Limit		Value	%	25		<table border="1"> <thead> <tr> <th colspan="2">Prozone Limit</th> </tr> <tr> <th>Value</th> <th>%</th> </tr> </thead> <tbody> <tr> <td>2500.0</td> <td></td> </tr> </tbody> </table>		Prozone Limit		Value	%	2500.0																																																																																
Non-Linear Limit																																																																																														
Value	%																																																																																													
25																																																																																														
Prozone Limit																																																																																														
Value	%																																																																																													
2500.0																																																																																														
<table border="1"> <thead> <tr> <th colspan="2">User Range</th> </tr> <tr> <th>Lower</th> <th>Upper</th> </tr> </thead> <tbody> <tr> <td>1</td> <td>1100</td> </tr> </tbody> </table>		User Range		Lower	Upper	1	1100	<table border="1"> <thead> <tr> <th colspan="2">Serum Index Limits</th> </tr> <tr> <th></th> <th>Limit</th> </tr> </thead> <tbody> <tr> <td>Hemolysis</td> <td>0.0</td> </tr> <tr> <td>Icterus</td> <td>0.0</td> </tr> <tr> <td>Lipemia</td> <td>0.0</td> </tr> </tbody> </table>		Serum Index Limits			Limit	Hemolysis	0.0	Icterus	0.0	Lipemia	0.0																																																																											
User Range																																																																																														
Lower	Upper																																																																																													
1	1100																																																																																													
Serum Index Limits																																																																																														
	Limit																																																																																													
Hemolysis	0.0																																																																																													
Icterus	0.0																																																																																													
Lipemia	0.0																																																																																													
<table border="1"> <thead> <tr> <th colspan="2">ValidRange</th> </tr> <tr> <th>Lower</th> <th>Upper</th> </tr> </thead> <tbody> <tr> <td>1</td> <td>1100</td> </tr> </tbody> </table>		ValidRange		Lower	Upper	1	1100	<table border="1"> <thead> <tr> <th colspan="2">Equation</th> </tr> <tr> <th>Limit</th> <th>Equation</th> </tr> </thead> <tbody> <tr> <td>0.0</td> <td>none</td> </tr> </tbody> </table>		Equation		Limit	Equation	0.0	none																																																																															
ValidRange																																																																																														
Lower	Upper																																																																																													
1	1100																																																																																													
Equation																																																																																														
Limit	Equation																																																																																													
0.0	none																																																																																													
<table border="1"> <thead> <tr> <th colspan="2">Qualitative</th> </tr> <tr> <th>On</th> <th>Off</th> </tr> </thead> <tbody> <tr> <td><input type="radio"/></td> <td><input checked="" type="radio"/></td> </tr> </tbody> </table>		Qualitative		On	Off	<input type="radio"/>	<input checked="" type="radio"/>	<table border="1"> <thead> <tr> <th colspan="2">Judge Point</th> </tr> <tr> <th>Value</th> <th>%</th> </tr> </thead> <tbody> <tr> <td>0</td> <td></td> </tr> </tbody> </table>		Judge Point		Value	%	0																																																																																
Qualitative																																																																																														
On	Off																																																																																													
<input type="radio"/>	<input checked="" type="radio"/>																																																																																													
Judge Point																																																																																														
Value	%																																																																																													
0																																																																																														

Test Name		Test Code		Measure		Sample Reagent		Ranges Limits		Calibration																																				
NN	CK	IL	CK-d																																											
Calibration				Condition				Calibration/R-Blank Limit																																						
<table border="1"> <thead> <tr> <th colspan="2">Method</th> <th colspan="2">Curve Type</th> </tr> </thead> <tbody> <tr> <td>1-point</td> <td></td> <td>Linear</td> <td></td> </tr> </tbody> </table>				Method		Curve Type		1-point		Linear		<table border="1"> <thead> <tr> <th colspan="2">Repeats</th> <th colspan="2">Stability</th> </tr> </thead> <tbody> <tr> <td>2</td> <td></td> <td></td> <td>days</td> </tr> </tbody> </table>				Repeats		Stability		2			days	<table border="1"> <thead> <tr> <th colspan="2">No</th> </tr> </thead> <tbody> <tr> <td></td> <td></td> </tr> </tbody> </table>				No																		
Method		Curve Type																																												
1-point		Linear																																												
Repeats		Stability																																												
2			days																																											
No																																														
<table border="1"> <thead> <tr> <th colspan="3">Calibrator</th> </tr> <tr> <th></th> <th>Calibrator</th> <th>Conc</th> </tr> </thead> <tbody> <tr> <td>1</td> <td>TruCalU</td> <td>**</td> </tr> <tr> <td>2</td> <td></td> <td></td> </tr> <tr> <td>3</td> <td></td> <td></td> </tr> <tr> <td>4</td> <td></td> <td></td> </tr> <tr> <td>5</td> <td></td> <td></td> </tr> </tbody> </table>				Calibrator				Calibrator	Conc	1	TruCalU	**	2			3			4			5			<table border="1"> <thead> <tr> <th colspan="2">Condition</th> </tr> </thead> <tbody> <tr> <td><input checked="" type="checkbox"/></td> <td>Reagent Blank</td> </tr> <tr> <td><input type="checkbox"/></td> <td>Auto Reagent Blank by Bottle</td> </tr> <tr> <td><input type="checkbox"/></td> <td>Ask for calibration when reagent lot changes</td> </tr> </tbody> </table>				Condition		<input checked="" type="checkbox"/>	Reagent Blank	<input type="checkbox"/>	Auto Reagent Blank by Bottle	<input type="checkbox"/>	Ask for calibration when reagent lot changes	<table border="1"> <thead> <tr> <th colspan="2">R-Blank Limit</th> </tr> <tr> <th>Value</th> <th>mAbs</th> </tr> </thead> <tbody> <tr> <td>2500.0</td> <td></td> </tr> </tbody> </table>				R-Blank Limit		Value	mAbs	2500.0	
Calibrator																																														
	Calibrator	Conc																																												
1	TruCalU	**																																												
2																																														
3																																														
4																																														
5																																														
Condition																																														
<input checked="" type="checkbox"/>	Reagent Blank																																													
<input type="checkbox"/>	Auto Reagent Blank by Bottle																																													
<input type="checkbox"/>	Ask for calibration when reagent lot changes																																													
R-Blank Limit																																														
Value	mAbs																																													
2500.0																																														
				<table border="1"> <thead> <tr> <th colspan="2">Cal Reps Range</th> </tr> <tr> <th>Value</th> <th>%</th> </tr> </thead> <tbody> <tr> <td>20.0</td> <td></td> </tr> </tbody> </table>				Cal Reps Range		Value	%	20.0																																		
Cal Reps Range																																														
Value	%																																													
20.0																																														
				<table border="1"> <thead> <tr> <th colspan="2">Min Cal Reps</th> </tr> <tr> <th>Value</th> <th>mAbs</th> </tr> </thead> <tbody> <tr> <td>50.0</td> <td></td> </tr> </tbody> </table>				Min Cal Reps		Value	mAbs	50.0																																		
Min Cal Reps																																														
Value	mAbs																																													
50.0																																														
				<table border="1"> <thead> <tr> <th colspan="2">Factor Change</th> </tr> <tr> <th>Value</th> <th>%</th> </tr> </thead> <tbody> <tr> <td>20</td> <td></td> </tr> </tbody> </table>				Factor Change		Value	%	20																																		
Factor Change																																														
Value	%																																													
20																																														
				<table border="1"> <thead> <tr> <th colspan="2">M-Point Curve Fit</th> </tr> <tr> <th>Value</th> <th>%</th> </tr> </thead> <tbody> <tr> <td>0.0</td> <td></td> </tr> </tbody> </table>				M-Point Curve Fit		Value	%	0.0																																		
M-Point Curve Fit																																														
Value	%																																													
0.0																																														

# Креатинкиназа

IFCC, DGKCOpt



\* -Рекомендуемый параметр.

**\*\*** - Вводится из паспорта калибратору.

Тип калибровки: Линейная (калибратор TruCalU) или по фактору (см инструкцию).

Контроль по TruLab N и TruLab P.