

Параметры для ввода в программу анализатора iLab-Taurus

Test Name: NN CK-MB Test Code: IL_ CK-MBd

Measure Sample Reagent Ranges Limits Calibration

Sample Type: Serum

Reaction Cycle: Standard Extended

Reporting Unit: Ед/л * Decimal Points: 1

User Define: Conv.Factor: 0.000

Methodology: Type: End Point Rate

Measuring Point: 27 - 33

Photometric: 1 Wavelength 2 Wavelength

Primary: 340 Secondary: 546

Correction Constant: Slope: 1.000 Intercept: 0.000

Test Name: NN CK-MB Test Code: IL_ CK-MBd

Measure Sample Reagent Ranges Limits Calibration

Sample Volume: Sampling Condition

Sample	Dilution	Diluent
Volume	Sample Vol.	Diluent Vol.
1 7.0	0.0	0.0
2 3.0	0.0	0.0
3 10.0	0.0	0.0
4 2.0	20	80

Diluent: Diluent Warning Limit: 0 Tests

Reagent Volume: R1 CK-MB1 R2 CK-MB2

	Volume	Diluent Vol.	stirring	Warning Limit(tests)	Stability(days)
R1	140	0.0	✓	20	No Control
R2	35	0.0	✓	20	No Control

	Condition No.	1	2	3	4
First Run		✓			
Samp.Vol.Reduction			✓		
Below N-Range					
Above N-Range					
Panic L					
Panic H			✓		
User Range L					
User Range H	✓		✓		
Noise					
Prozone					
HIGH!	✓		✓		
ABS!	✓		✓		

Ranges		Limits																																																																																												
<table border="1"> <thead> <tr> <th colspan="3">Normal Range</th> </tr> <tr> <th></th> <th>Lower</th> <th>Upper</th> </tr> </thead> <tbody> <tr> <td>Male</td> <td>0</td> <td>24</td> </tr> <tr> <td>Female</td> <td>0</td> <td>24</td> </tr> <tr> <td>Other</td> <td>0</td> <td>24</td> </tr> </tbody> </table>		Normal Range				Lower	Upper	Male	0	24	Female	0	24	Other	0	24	<table border="1"> <thead> <tr> <th colspan="2">Reaction Slope</th> <th colspan="2">Absorbance Limit</th> </tr> </thead> <tbody> <tr> <td colspan="2"> <input type="radio"/> Negative <input checked="" type="radio"/> Positive </td> <td colspan="2"> <input checked="" type="radio"/> Above <input type="radio"/> Below </td> </tr> <tr> <td colspan="2"> <table border="1"> <thead> <tr> <th colspan="2">Non-Linear Limit</th> </tr> <tr> <th>Value</th> <th>%</th> </tr> </thead> <tbody> <tr> <td>25</td> <td></td> </tr> </tbody> </table> </td> <td colspan="2"> <table border="1"> <thead> <tr> <th colspan="2">Prozone Limit</th> </tr> <tr> <th>Value</th> <th>%</th> </tr> </thead> <tbody> <tr> <td>2500.0</td> <td></td> </tr> </tbody> </table> </td> </tr> <tr> <td colspan="2"> <table border="1"> <thead> <tr> <th colspan="2">User Range</th> </tr> <tr> <th>Lower</th> <th>Upper</th> </tr> </thead> <tbody> <tr> <td>1</td> <td>2000</td> </tr> </tbody> </table> </td> <td colspan="2"> <table border="1"> <thead> <tr> <th colspan="2">Serum Index Limits</th> </tr> <tr> <th></th> <th>Limit</th> </tr> </thead> <tbody> <tr> <td>Hemolysis</td> <td>0.0</td> </tr> <tr> <td>Icterus</td> <td>0.0</td> </tr> <tr> <td>Lipemia</td> <td>0.0</td> </tr> </tbody> </table> </td> </tr> <tr> <td colspan="2"> <table border="1"> <thead> <tr> <th colspan="2">ValidRange</th> </tr> <tr> <th>Lower</th> <th>Upper</th> </tr> </thead> <tbody> <tr> <td>1</td> <td>2000</td> </tr> </tbody> </table> </td> <td colspan="2"> <table border="1"> <thead> <tr> <th colspan="2">Equation</th> </tr> <tr> <th>Limit</th> <th>Equation</th> </tr> </thead> <tbody> <tr> <td>0.0</td> <td>none</td> </tr> </tbody> </table> </td> </tr> <tr> <td colspan="2"> <table border="1"> <thead> <tr> <th colspan="2">Qualitative</th> </tr> <tr> <th>On</th> <th>Off</th> </tr> </thead> <tbody> <tr> <td><input type="radio"/></td> <td><input checked="" type="radio"/></td> </tr> </tbody> </table> </td> <td colspan="2"> <table border="1"> <thead> <tr> <th colspan="2">Judge Point</th> </tr> <tr> <th>Value</th> <th>%</th> </tr> </thead> <tbody> <tr> <td>0</td> <td></td> </tr> </tbody> </table> </td> </tr> </tbody> </table>		Reaction Slope		Absorbance Limit		<input type="radio"/> Negative <input checked="" type="radio"/> Positive		<input checked="" type="radio"/> Above <input type="radio"/> Below		<table border="1"> <thead> <tr> <th colspan="2">Non-Linear Limit</th> </tr> <tr> <th>Value</th> <th>%</th> </tr> </thead> <tbody> <tr> <td>25</td> <td></td> </tr> </tbody> </table>		Non-Linear Limit		Value	%	25		<table border="1"> <thead> <tr> <th colspan="2">Prozone Limit</th> </tr> <tr> <th>Value</th> <th>%</th> </tr> </thead> <tbody> <tr> <td>2500.0</td> <td></td> </tr> </tbody> </table>		Prozone Limit		Value	%	2500.0		<table border="1"> <thead> <tr> <th colspan="2">User Range</th> </tr> <tr> <th>Lower</th> <th>Upper</th> </tr> </thead> <tbody> <tr> <td>1</td> <td>2000</td> </tr> </tbody> </table>		User Range		Lower	Upper	1	2000	<table border="1"> <thead> <tr> <th colspan="2">Serum Index Limits</th> </tr> <tr> <th></th> <th>Limit</th> </tr> </thead> <tbody> <tr> <td>Hemolysis</td> <td>0.0</td> </tr> <tr> <td>Icterus</td> <td>0.0</td> </tr> <tr> <td>Lipemia</td> <td>0.0</td> </tr> </tbody> </table>		Serum Index Limits			Limit	Hemolysis	0.0	Icterus	0.0	Lipemia	0.0	<table border="1"> <thead> <tr> <th colspan="2">ValidRange</th> </tr> <tr> <th>Lower</th> <th>Upper</th> </tr> </thead> <tbody> <tr> <td>1</td> <td>2000</td> </tr> </tbody> </table>		ValidRange		Lower	Upper	1	2000	<table border="1"> <thead> <tr> <th colspan="2">Equation</th> </tr> <tr> <th>Limit</th> <th>Equation</th> </tr> </thead> <tbody> <tr> <td>0.0</td> <td>none</td> </tr> </tbody> </table>		Equation		Limit	Equation	0.0	none	<table border="1"> <thead> <tr> <th colspan="2">Qualitative</th> </tr> <tr> <th>On</th> <th>Off</th> </tr> </thead> <tbody> <tr> <td><input type="radio"/></td> <td><input checked="" type="radio"/></td> </tr> </tbody> </table>		Qualitative		On	Off	<input type="radio"/>	<input checked="" type="radio"/>	<table border="1"> <thead> <tr> <th colspan="2">Judge Point</th> </tr> <tr> <th>Value</th> <th>%</th> </tr> </thead> <tbody> <tr> <td>0</td> <td></td> </tr> </tbody> </table>		Judge Point		Value	%	0	
Normal Range																																																																																														
	Lower	Upper																																																																																												
Male	0	24																																																																																												
Female	0	24																																																																																												
Other	0	24																																																																																												
Reaction Slope		Absorbance Limit																																																																																												
<input type="radio"/> Negative <input checked="" type="radio"/> Positive		<input checked="" type="radio"/> Above <input type="radio"/> Below																																																																																												
<table border="1"> <thead> <tr> <th colspan="2">Non-Linear Limit</th> </tr> <tr> <th>Value</th> <th>%</th> </tr> </thead> <tbody> <tr> <td>25</td> <td></td> </tr> </tbody> </table>		Non-Linear Limit		Value	%	25		<table border="1"> <thead> <tr> <th colspan="2">Prozone Limit</th> </tr> <tr> <th>Value</th> <th>%</th> </tr> </thead> <tbody> <tr> <td>2500.0</td> <td></td> </tr> </tbody> </table>		Prozone Limit		Value	%	2500.0																																																																																
Non-Linear Limit																																																																																														
Value	%																																																																																													
25																																																																																														
Prozone Limit																																																																																														
Value	%																																																																																													
2500.0																																																																																														
<table border="1"> <thead> <tr> <th colspan="2">User Range</th> </tr> <tr> <th>Lower</th> <th>Upper</th> </tr> </thead> <tbody> <tr> <td>1</td> <td>2000</td> </tr> </tbody> </table>		User Range		Lower	Upper	1	2000	<table border="1"> <thead> <tr> <th colspan="2">Serum Index Limits</th> </tr> <tr> <th></th> <th>Limit</th> </tr> </thead> <tbody> <tr> <td>Hemolysis</td> <td>0.0</td> </tr> <tr> <td>Icterus</td> <td>0.0</td> </tr> <tr> <td>Lipemia</td> <td>0.0</td> </tr> </tbody> </table>		Serum Index Limits			Limit	Hemolysis	0.0	Icterus	0.0	Lipemia	0.0																																																																											
User Range																																																																																														
Lower	Upper																																																																																													
1	2000																																																																																													
Serum Index Limits																																																																																														
	Limit																																																																																													
Hemolysis	0.0																																																																																													
Icterus	0.0																																																																																													
Lipemia	0.0																																																																																													
<table border="1"> <thead> <tr> <th colspan="2">ValidRange</th> </tr> <tr> <th>Lower</th> <th>Upper</th> </tr> </thead> <tbody> <tr> <td>1</td> <td>2000</td> </tr> </tbody> </table>		ValidRange		Lower	Upper	1	2000	<table border="1"> <thead> <tr> <th colspan="2">Equation</th> </tr> <tr> <th>Limit</th> <th>Equation</th> </tr> </thead> <tbody> <tr> <td>0.0</td> <td>none</td> </tr> </tbody> </table>		Equation		Limit	Equation	0.0	none																																																																															
ValidRange																																																																																														
Lower	Upper																																																																																													
1	2000																																																																																													
Equation																																																																																														
Limit	Equation																																																																																													
0.0	none																																																																																													
<table border="1"> <thead> <tr> <th colspan="2">Qualitative</th> </tr> <tr> <th>On</th> <th>Off</th> </tr> </thead> <tbody> <tr> <td><input type="radio"/></td> <td><input checked="" type="radio"/></td> </tr> </tbody> </table>		Qualitative		On	Off	<input type="radio"/>	<input checked="" type="radio"/>	<table border="1"> <thead> <tr> <th colspan="2">Judge Point</th> </tr> <tr> <th>Value</th> <th>%</th> </tr> </thead> <tbody> <tr> <td>0</td> <td></td> </tr> </tbody> </table>		Judge Point		Value	%	0																																																																																
Qualitative																																																																																														
On	Off																																																																																													
<input type="radio"/>	<input checked="" type="radio"/>																																																																																													
Judge Point																																																																																														
Value	%																																																																																													
0																																																																																														

Test Name		Test Code		Measure		Sample Reagent		Ranges Limits		Calibration													
NN	CK-MB	IL	CK-MBd																				
Calibration			Condition			Calibration/R-Blank Limit																	
<table border="1"> <thead> <tr> <th colspan="2">Method</th> </tr> <tr> <th>Value</th> <th>Unit</th> </tr> </thead> <tbody> <tr> <td>Factor</td> <td></td> </tr> </tbody> </table>			Method		Value	Unit	Factor		<input checked="" type="checkbox"/> Reagent Blank <input type="checkbox"/> Auto Reagent Blank by Bottle <input type="checkbox"/> Ask for calibration when reagent lot changes			<table border="1"> <thead> <tr> <th colspan="2">R-Blank Limit</th> </tr> <tr> <th>Value</th> <th>Unit</th> </tr> </thead> <tbody> <tr> <td>2500.0</td> <td>mAbs</td> </tr> </tbody> </table>						R-Blank Limit		Value	Unit	2500.0	mAbs
Method																							
Value	Unit																						
Factor																							
R-Blank Limit																							
Value	Unit																						
2500.0	mAbs																						
						<table border="1"> <thead> <tr> <th colspan="2">Cal Reps Range</th> </tr> <tr> <th>Value</th> <th>%</th> </tr> </thead> <tbody> <tr> <td></td> <td></td> </tr> </tbody> </table>						Cal Reps Range		Value	%								
Cal Reps Range																							
Value	%																						
						<table border="1"> <thead> <tr> <th colspan="2">Min Cal Reps</th> </tr> <tr> <th>Value</th> <th>Unit</th> </tr> </thead> <tbody> <tr> <td></td> <td>mAbs</td> </tr> </tbody> </table>						Min Cal Reps		Value	Unit		mAbs						
Min Cal Reps																							
Value	Unit																						
	mAbs																						
						<table border="1"> <thead> <tr> <th colspan="2">Factor Change</th> </tr> <tr> <th>Value</th> <th>%</th> </tr> </thead> <tbody> <tr> <td>0</td> <td></td> </tr> </tbody> </table>						Factor Change		Value	%	0							
Factor Change																							
Value	%																						
0																							
						<table border="1"> <thead> <tr> <th colspan="2">M-Point Curve Fit</th> </tr> <tr> <th>Value</th> <th>%</th> </tr> </thead> <tbody> <tr> <td>0.0</td> <td></td> </tr> </tbody> </table>						M-Point Curve Fit		Value	%	0.0							
M-Point Curve Fit																							
Value	%																						
0.0																							
						<table border="1"> <thead> <tr> <th colspan="2">Factor</th> </tr> <tr> <th>Value</th> <th>Unit</th> </tr> </thead> <tbody> <tr> <td>9000</td> <td></td> </tr> </tbody> </table>						Factor		Value	Unit	9000							
Factor																							
Value	Unit																						
9000																							

* -Рекомендуемый параметр.

Тип калибровки: по фактору (см инструкцию). Можно использовать сторонний калибратор (Например фирмы Roche), в таком случае калибровка 1-точечная линейная.

Контроль по TruLab N и TruLab P.