

Параметры для ввода в программу анализатора Beckman-Coulter AU 480

Parameters		Specific Test Parameters																																																																																																																																																																																																									
General	LIH	ISE	HbA1c	Calculated Test	Range																																																																																																																																																																																																						
Test Name: <input type="text" value="ASLO"/> < > Type: <input type="text" value="Serum"/> Operation <input type="text" value="Yes"/>																																																																																																																																																																																																											
<table style="width: 100%; border: none;"> <tr> <td style="width: 20%;">Sample Volume</td> <td style="width: 10%;"><input type="text" value="1.2"/></td> <td style="width: 10%;">μL</td> <td style="width: 10%;">Dilution</td> <td style="width: 10%;"><input type="text" value="0"/></td> <td style="width: 10%;">μL</td> <td style="width: 10%;">OD Limit</td> <td colspan="2"></td> </tr> <tr> <td>Pre-Dilution Rate</td> <td><input type="text" value="1"/></td> <td></td> <td></td> <td></td> <td></td> <td>Min.OD</td> <td><input type="text" value="-2.0"/></td> <td>Max.OD</td> <td><input type="text" value="2.5"/></td> </tr> <tr> <td>Rgt. Volume</td> <td>R1(R1-1)</td> <td><input type="text" value="100"/></td> <td>μL</td> <td>Dilution</td> <td><input type="text" value="0"/></td> <td>μL</td> <td>Reagent OD Limit</td> <td colspan="2"></td> </tr> <tr> <td></td> <td></td> <td></td> <td></td> <td></td> <td></td> <td></td> <td>First</td> <td>Low</td> <td><input type="text" value="-2.0"/></td> <td>High</td> <td><input type="text" value="2.5"/></td> </tr> <tr> <td></td> <td></td> <td></td> <td></td> <td></td> <td></td> <td></td> <td>Last</td> <td>Low</td> <td><input type="text" value="-2.0"/></td> <td>High</td> <td><input type="text" value="2.5"/></td> </tr> <tr> <td></td> <td>R2(R2-1)</td> <td><input type="text" value="20"/></td> <td>μL</td> <td>Dilution</td> <td><input type="text" value="0"/></td> <td>μL</td> <td>Dynamic Range Low</td> <td><input type="text" value="5"/></td> <td>High</td> <td><input type="text" value="800"/></td> </tr> <tr> <td>Wavelength</td> <td>Pri</td> <td><input type="text" value="540"/></td> <td>nm</td> <td>Sec.</td> <td><input type="text" value="NONE"/></td> <td>nm</td> <td>Correlation Factor A</td> <td><input type="text" value="1"/></td> <td>B</td> <td><input type="text" value="0"/></td> </tr> <tr> <td>Method</td> <td></td> <td><input type="text" value="RATE"/></td> <td></td> <td></td> <td></td> <td></td> <td>Factor for Maker</td> <td>A</td> <td><input type="text" value="1"/></td> <td>B</td> <td><input type="text" value="0"/></td> </tr> <tr> <td>Reaction Slope</td> <td></td> <td><input type="text" value="+"/></td> <td></td> <td></td> <td></td> <td></td> <td>Onboard Stability Period</td> <td><input type="text" value="999"/></td> <td>Day</td> <td><input type="text" value=""/></td> <td>Hour</td> </tr> <tr> <td>Measuring Point1</td> <td>First</td> <td><input type="text" value="13"/></td> <td></td> <td>Last</td> <td><input type="text" value="23"/></td> <td></td> <td colspan="4"></td> </tr> <tr> <td>Measuring Point2</td> <td>First</td> <td><input type="text" value=""/></td> <td></td> <td>Last</td> <td><input type="text" value=""/></td> <td></td> <td colspan="4"></td> </tr> <tr> <td>Linearity Limit</td> <td></td> <td><input type="text" value=""/></td> <td>%</td> <td colspan="7"></td> </tr> <tr> <td>Lag Time Check</td> <td></td> <td><input type="text" value=""/></td> <td></td> <td colspan="7"></td> </tr> </table>						Sample Volume	<input type="text" value="1.2"/>	μL	Dilution	<input type="text" value="0"/>	μL	OD Limit			Pre-Dilution Rate	<input type="text" value="1"/>					Min.OD	<input type="text" value="-2.0"/>	Max.OD	<input type="text" value="2.5"/>	Rgt. Volume	R1(R1-1)	<input type="text" value="100"/>	μL	Dilution	<input type="text" value="0"/>	μL	Reagent OD Limit										First	Low	<input type="text" value="-2.0"/>	High	<input type="text" value="2.5"/>								Last	Low	<input type="text" value="-2.0"/>	High	<input type="text" value="2.5"/>		R2(R2-1)	<input type="text" value="20"/>	μL	Dilution	<input type="text" value="0"/>	μL	Dynamic Range Low	<input type="text" value="5"/>	High	<input type="text" value="800"/>	Wavelength	Pri	<input type="text" value="540"/>	nm	Sec.	<input type="text" value="NONE"/>	nm	Correlation Factor A	<input type="text" value="1"/>	B	<input type="text" value="0"/>	Method		<input type="text" value="RATE"/>					Factor for Maker	A	<input type="text" value="1"/>	B	<input type="text" value="0"/>	Reaction Slope		<input type="text" value="+"/>					Onboard Stability Period	<input type="text" value="999"/>	Day	<input type="text" value=""/>	Hour	Measuring Point1	First	<input type="text" value="13"/>		Last	<input type="text" value="23"/>						Measuring Point2	First	<input type="text" value=""/>		Last	<input type="text" value=""/>						Linearity Limit		<input type="text" value=""/>	%								Lag Time Check		<input type="text" value=""/>																																																															
Sample Volume	<input type="text" value="1.2"/>	μL	Dilution	<input type="text" value="0"/>	μL	OD Limit																																																																																																																																																																																																					
Pre-Dilution Rate	<input type="text" value="1"/>					Min.OD	<input type="text" value="-2.0"/>	Max.OD	<input type="text" value="2.5"/>																																																																																																																																																																																																		
Rgt. Volume	R1(R1-1)	<input type="text" value="100"/>	μL	Dilution	<input type="text" value="0"/>	μL	Reagent OD Limit																																																																																																																																																																																																				
							First	Low	<input type="text" value="-2.0"/>	High	<input type="text" value="2.5"/>																																																																																																																																																																																																
							Last	Low	<input type="text" value="-2.0"/>	High	<input type="text" value="2.5"/>																																																																																																																																																																																																
	R2(R2-1)	<input type="text" value="20"/>	μL	Dilution	<input type="text" value="0"/>	μL	Dynamic Range Low	<input type="text" value="5"/>	High	<input type="text" value="800"/>																																																																																																																																																																																																	
Wavelength	Pri	<input type="text" value="540"/>	nm	Sec.	<input type="text" value="NONE"/>	nm	Correlation Factor A	<input type="text" value="1"/>	B	<input type="text" value="0"/>																																																																																																																																																																																																	
Method		<input type="text" value="RATE"/>					Factor for Maker	A	<input type="text" value="1"/>	B	<input type="text" value="0"/>																																																																																																																																																																																																
Reaction Slope		<input type="text" value="+"/>					Onboard Stability Period	<input type="text" value="999"/>	Day	<input type="text" value=""/>	Hour																																																																																																																																																																																																
Measuring Point1	First	<input type="text" value="13"/>		Last	<input type="text" value="23"/>																																																																																																																																																																																																						
Measuring Point2	First	<input type="text" value=""/>		Last	<input type="text" value=""/>																																																																																																																																																																																																						
Linearity Limit		<input type="text" value=""/>	%																																																																																																																																																																																																								
Lag Time Check		<input type="text" value=""/>																																																																																																																																																																																																									
<table style="width: 100%; border: none;"> <tr> <td style="width: 15%;">Value/Flag:</td> <td style="width: 10%;"><input type="text" value=""/></td> <td style="width: 10%;">Level L:</td> <td style="width: 10%;"><input type="text" value=""/></td> <td style="width: 10%;">Level H:</td> <td style="width: 10%;"><input type="text" value=""/></td> <td colspan="5"></td> </tr> <tr> <td colspan="6">Specific Ranges:</td> <td colspan="5" style="border: 1px solid black; padding: 5px;">Panic Value</td> </tr> <tr> <td></td> <td>Sex</td> <td>Year</td> <td>Month</td> <td>Year</td> <td>Month</td> <td>Low</td> <td>High</td> <td colspan="3"></td> </tr> <tr> <td><input type="checkbox"/></td> <td>1.</td> <td><input type="text" value=""/></td> <td colspan="3"></td> </tr> <tr> <td><input type="checkbox"/></td> <td>2.</td> <td><input type="text" value=""/></td> <td colspan="3"></td> </tr> <tr> <td><input type="checkbox"/></td> <td>3.</td> <td><input type="text" value=""/></td> <td colspan="3"></td> </tr> <tr> <td><input type="checkbox"/></td> <td>4.</td> <td><input type="text" value=""/></td> <td colspan="3"></td> </tr> <tr> <td><input type="checkbox"/></td> <td>5.</td> <td><input type="text" value=""/></td> <td colspan="3"></td> </tr> <tr> <td><input type="checkbox"/></td> <td>6.</td> <td><input type="text" value=""/></td> <td colspan="3"></td> </tr> <tr> <td><input type="checkbox"/></td> <td colspan="5">7. No demographics</td> <td><input type="text" value="0"/></td> <td><input type="text" value="200"/></td> <td colspan="3"></td> </tr> <tr> <td><input type="checkbox"/></td> <td colspan="5">8. Not within expected values</td> <td><input type="text" value="0"/></td> <td><input type="text" value="200"/></td> <td colspan="3"></td> </tr> <tr> <td>Unit</td> <td colspan="5"><input type="text" value="ME/ml"/></td> <td>Decimal Places</td> <td><input type="text" value="0"/></td> <td colspan="3"></td> </tr> </table>						Value/Flag:	<input type="text" value=""/>	Level L:	<input type="text" value=""/>	Level H:	<input type="text" value=""/>						Specific Ranges:						Panic Value						Sex	Year	Month	Year	Month	Low	High				<input type="checkbox"/>	1.	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>				<input type="checkbox"/>	2.	<input type="text" value=""/>	<input type="text" value=""/>				<input type="checkbox"/>	3.	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>				<input type="checkbox"/>	4.	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>				<input type="checkbox"/>	5.	<input type="text" value=""/>				<input type="checkbox"/>	6.	<input type="text" value=""/>				<input type="checkbox"/>	7. No demographics					<input type="text" value="0"/>	<input type="text" value="200"/>				<input type="checkbox"/>	8. Not within expected values					<input type="text" value="0"/>	<input type="text" value="200"/>				Unit	<input type="text" value="ME/ml"/>					Decimal Places	<input type="text" value="0"/>																																																																																			
Value/Flag:	<input type="text" value=""/>	Level L:	<input type="text" value=""/>	Level H:	<input type="text" value=""/>																																																																																																																																																																																																						
Specific Ranges:						Panic Value																																																																																																																																																																																																					
	Sex	Year	Month	Year	Month	Low	High																																																																																																																																																																																																				
<input type="checkbox"/>	1.	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>																																																																																																																																																																																																				
<input type="checkbox"/>	2.	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>																																																																																																																																																																																																				
<input type="checkbox"/>	3.	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>																																																																																																																																																																																																				
<input type="checkbox"/>	4.	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>																																																																																																																																																																																																				
<input type="checkbox"/>	5.	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>																																																																																																																																																																																																				
<input type="checkbox"/>	6.	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>																																																																																																																																																																																																				
<input type="checkbox"/>	7. No demographics					<input type="text" value="0"/>	<input type="text" value="200"/>																																																																																																																																																																																																				
<input type="checkbox"/>	8. Not within expected values					<input type="text" value="0"/>	<input type="text" value="200"/>																																																																																																																																																																																																				
Unit	<input type="text" value="ME/ml"/>					Decimal Places	<input type="text" value="0"/>																																																																																																																																																																																																				
<table style="width: 100%; border: none;"> <tr> <th colspan="2" style="text-align: left;">Parameters</th> <th colspan="4" style="text-align: right;">Calibration Parameters</th> </tr> <tr> <th style="width: 15%;">Calibrators</th> <th style="width: 10%;">Calibration Specific</th> <th colspan="4" style="text-align: right;">STAT Table Calibration</th> </tr> <tr> <td colspan="6">General ISE</td> </tr> <tr> <td colspan="6"> Test Name: <input type="text" value="ASLO"/> < > Type: <input type="text" value="Serum"/> <input type="checkbox"/> Use Serum Cal. </td> </tr> <tr> <td colspan="6"> <table style="width: 100%; border: none;"> <tr> <td style="width: 20%;">Calibration Type:</td> <td style="width: 10%;"><input type="text" value="6AB"/></td> <td style="width: 10%;">Formula:</td> <td style="width: 10%;"><input type="text" value="Spline"/></td> <td style="width: 10%;">Counts:</td> <td style="width: 10%;"><input type="text" value="2"/></td> <td colspan="5"></td> </tr> <tr> <td colspan="6"><Calibrator Parameters></td> </tr> <tr> <td>Calibrator</td> <td>OD</td> <td>Conc</td> <td>Factor/OD range</td> <td>Low</td> <td>High</td> <td colspan="5">Slope Check <input type="text" value="None"/></td> </tr> <tr> <td>Point 1:</td> <td>Saline</td> <td></td> <td></td> <td>-0.1</td> <td>2.5</td> <td colspan="5" rowspan="5"> Allowance Range Check <input type="checkbox"/> Reagent Blank <input type="checkbox"/> Calibration Advanced Calibration Operation <input type="text" value=""/> Interval (RB/ACAL) <input type="text" value=""/> </td> </tr> <tr> <td>Point 2:</td> <td>TruCal ASLO 1</td> <td></td> <td>*</td> <td>-0.1</td> <td>2.5</td> </tr> <tr> <td>Point 3:</td> <td>TruCal ASLO 2</td> <td></td> <td>*</td> <td>-0.1</td> <td>2.5</td> </tr> <tr> <td>Point 4:</td> <td>TruCal ASLO 3</td> <td></td> <td>*</td> <td>-0.1</td> <td>2.5</td> </tr> <tr> <td>Point 5:</td> <td>TruCal ASLO 4</td> <td></td> <td>*</td> <td>-0.1</td> <td>2.5</td> </tr> <tr> <td>Point 6:</td> <td>TruCal ASLO 5</td> <td></td> <td>*</td> <td>-0.1</td> <td>2.5</td> <td colspan="5"></td> </tr> <tr> <td>Point 7:</td> <td></td> <td></td> <td></td> <td></td> <td></td> <td colspan="5"></td> </tr> <tr> <td>Point 8:</td> <td></td> <td></td> <td></td> <td></td> <td></td> <td colspan="5"></td> </tr> <tr> <td>Point 9:</td> <td></td> <td></td> <td></td> <td></td> <td></td> <td colspan="5"></td> </tr> <tr> <td>Point 10:</td> <td></td> <td></td> <td></td> <td></td> <td></td> <td colspan="5"></td> </tr> <tr> <td colspan="6"><Point Cal. For Master Curve></td> </tr> <tr> <td>Calibrator</td> <td>OD</td> <td>Conc</td> <td>Low</td> <td>High</td> <td>Stability</td> <td colspan="5"></td> </tr> <tr> <td>Point 1:</td> <td></td> <td></td> <td></td> <td></td> <td>Reagent Blank</td> <td><input type="text" value=""/></td> <td>Day</td> <td><input type="text" value=""/></td> <td>Hour</td> <td></td> </tr> <tr> <td>Point 2:</td> <td></td> <td></td> <td></td> <td></td> <td>Calibration</td> <td><input type="text" value=""/></td> <td>Day</td> <td><input type="text" value=""/></td> <td>Hour</td> <td></td> </tr> <tr> <td>MB Type Factor:</td> <td colspan="5"><input type="text" value=""/></td> <td>1-Point Calibration Point</td> <td><input type="text" value=""/></td> <td colspan="3"><input type="checkbox"/> with Conc-0</td> </tr> </table> </td> </tr> </table>						Parameters		Calibration Parameters				Calibrators	Calibration Specific	STAT Table Calibration				General ISE						Test Name: <input type="text" value="ASLO"/> < > Type: <input type="text" value="Serum"/> <input type="checkbox"/> Use Serum Cal.						<table style="width: 100%; border: none;"> <tr> <td style="width: 20%;">Calibration Type:</td> <td style="width: 10%;"><input type="text" value="6AB"/></td> <td style="width: 10%;">Formula:</td> <td style="width: 10%;"><input type="text" value="Spline"/></td> <td style="width: 10%;">Counts:</td> <td style="width: 10%;"><input type="text" value="2"/></td> <td colspan="5"></td> </tr> <tr> <td colspan="6"><Calibrator Parameters></td> </tr> <tr> <td>Calibrator</td> <td>OD</td> <td>Conc</td> <td>Factor/OD range</td> <td>Low</td> <td>High</td> <td colspan="5">Slope Check <input type="text" value="None"/></td> </tr> <tr> <td>Point 1:</td> <td>Saline</td> <td></td> <td></td> <td>-0.1</td> <td>2.5</td> <td colspan="5" rowspan="5"> Allowance Range Check <input type="checkbox"/> Reagent Blank <input type="checkbox"/> Calibration Advanced Calibration Operation <input type="text" value=""/> Interval (RB/ACAL) <input type="text" value=""/> </td> </tr> <tr> <td>Point 2:</td> <td>TruCal ASLO 1</td> <td></td> <td>*</td> <td>-0.1</td> <td>2.5</td> </tr> <tr> <td>Point 3:</td> <td>TruCal ASLO 2</td> <td></td> <td>*</td> <td>-0.1</td> <td>2.5</td> </tr> <tr> <td>Point 4:</td> <td>TruCal ASLO 3</td> <td></td> <td>*</td> <td>-0.1</td> <td>2.5</td> </tr> <tr> <td>Point 5:</td> <td>TruCal ASLO 4</td> <td></td> <td>*</td> <td>-0.1</td> <td>2.5</td> </tr> <tr> <td>Point 6:</td> <td>TruCal ASLO 5</td> <td></td> <td>*</td> <td>-0.1</td> <td>2.5</td> <td colspan="5"></td> </tr> <tr> <td>Point 7:</td> <td></td> <td></td> <td></td> <td></td> <td></td> <td colspan="5"></td> </tr> <tr> <td>Point 8:</td> <td></td> <td></td> <td></td> <td></td> <td></td> <td colspan="5"></td> </tr> <tr> <td>Point 9:</td> <td></td> <td></td> <td></td> <td></td> <td></td> <td colspan="5"></td> </tr> <tr> <td>Point 10:</td> <td></td> <td></td> <td></td> <td></td> <td></td> <td colspan="5"></td> </tr> <tr> <td colspan="6"><Point Cal. For Master Curve></td> </tr> <tr> <td>Calibrator</td> <td>OD</td> <td>Conc</td> <td>Low</td> <td>High</td> <td>Stability</td> <td colspan="5"></td> </tr> <tr> <td>Point 1:</td> <td></td> <td></td> <td></td> <td></td> <td>Reagent Blank</td> <td><input type="text" value=""/></td> <td>Day</td> <td><input type="text" value=""/></td> <td>Hour</td> <td></td> </tr> <tr> <td>Point 2:</td> <td></td> <td></td> <td></td> <td></td> <td>Calibration</td> <td><input type="text" value=""/></td> <td>Day</td> <td><input type="text" value=""/></td> <td>Hour</td> <td></td> </tr> <tr> <td>MB Type Factor:</td> <td colspan="5"><input type="text" value=""/></td> <td>1-Point Calibration Point</td> <td><input type="text" value=""/></td> <td colspan="3"><input type="checkbox"/> with Conc-0</td> </tr> </table>						Calibration Type:	<input type="text" value="6AB"/>	Formula:	<input type="text" value="Spline"/>	Counts:	<input type="text" value="2"/>						<Calibrator Parameters>						Calibrator	OD	Conc	Factor/OD range	Low	High	Slope Check <input type="text" value="None"/>					Point 1:	Saline			-0.1	2.5	Allowance Range Check <input type="checkbox"/> Reagent Blank <input type="checkbox"/> Calibration Advanced Calibration Operation <input type="text" value=""/> Interval (RB/ACAL) <input type="text" value=""/>					Point 2:	TruCal ASLO 1		*	-0.1	2.5	Point 3:	TruCal ASLO 2		*	-0.1	2.5	Point 4:	TruCal ASLO 3		*	-0.1	2.5	Point 5:	TruCal ASLO 4		*	-0.1	2.5	Point 6:	TruCal ASLO 5		*	-0.1	2.5						Point 7:											Point 8:											Point 9:											Point 10:											<Point Cal. For Master Curve>						Calibrator	OD	Conc	Low	High	Stability						Point 1:					Reagent Blank	<input type="text" value=""/>	Day	<input type="text" value=""/>	Hour		Point 2:					Calibration	<input type="text" value=""/>	Day	<input type="text" value=""/>	Hour		MB Type Factor:	<input type="text" value=""/>					1-Point Calibration Point	<input type="text" value=""/>	<input type="checkbox"/> with Conc-0		
Parameters		Calibration Parameters																																																																																																																																																																																																									
Calibrators	Calibration Specific	STAT Table Calibration																																																																																																																																																																																																									
General ISE																																																																																																																																																																																																											
Test Name: <input type="text" value="ASLO"/> < > Type: <input type="text" value="Serum"/> <input type="checkbox"/> Use Serum Cal.																																																																																																																																																																																																											
<table style="width: 100%; border: none;"> <tr> <td style="width: 20%;">Calibration Type:</td> <td style="width: 10%;"><input type="text" value="6AB"/></td> <td style="width: 10%;">Formula:</td> <td style="width: 10%;"><input type="text" value="Spline"/></td> <td style="width: 10%;">Counts:</td> <td style="width: 10%;"><input type="text" value="2"/></td> <td colspan="5"></td> </tr> <tr> <td colspan="6"><Calibrator Parameters></td> </tr> <tr> <td>Calibrator</td> <td>OD</td> <td>Conc</td> <td>Factor/OD range</td> <td>Low</td> <td>High</td> <td colspan="5">Slope Check <input type="text" value="None"/></td> </tr> <tr> <td>Point 1:</td> <td>Saline</td> <td></td> <td></td> <td>-0.1</td> <td>2.5</td> <td colspan="5" rowspan="5"> Allowance Range Check <input type="checkbox"/> Reagent Blank <input type="checkbox"/> Calibration Advanced Calibration Operation <input type="text" value=""/> Interval (RB/ACAL) <input type="text" value=""/> </td> </tr> <tr> <td>Point 2:</td> <td>TruCal ASLO 1</td> <td></td> <td>*</td> <td>-0.1</td> <td>2.5</td> </tr> <tr> <td>Point 3:</td> <td>TruCal ASLO 2</td> <td></td> <td>*</td> <td>-0.1</td> <td>2.5</td> </tr> <tr> <td>Point 4:</td> <td>TruCal ASLO 3</td> <td></td> <td>*</td> <td>-0.1</td> <td>2.5</td> </tr> <tr> <td>Point 5:</td> <td>TruCal ASLO 4</td> <td></td> <td>*</td> <td>-0.1</td> <td>2.5</td> </tr> <tr> <td>Point 6:</td> <td>TruCal ASLO 5</td> <td></td> <td>*</td> <td>-0.1</td> <td>2.5</td> <td colspan="5"></td> </tr> <tr> <td>Point 7:</td> <td></td> <td></td> <td></td> <td></td> <td></td> <td colspan="5"></td> </tr> <tr> <td>Point 8:</td> <td></td> <td></td> <td></td> <td></td> <td></td> <td colspan="5"></td> </tr> <tr> <td>Point 9:</td> <td></td> <td></td> <td></td> <td></td> <td></td> <td colspan="5"></td> </tr> <tr> <td>Point 10:</td> <td></td> <td></td> <td></td> <td></td> <td></td> <td colspan="5"></td> </tr> <tr> <td colspan="6"><Point Cal. For Master Curve></td> </tr> <tr> <td>Calibrator</td> <td>OD</td> <td>Conc</td> <td>Low</td> <td>High</td> <td>Stability</td> <td colspan="5"></td> </tr> <tr> <td>Point 1:</td> <td></td> <td></td> <td></td> <td></td> <td>Reagent Blank</td> <td><input type="text" value=""/></td> <td>Day</td> <td><input type="text" value=""/></td> <td>Hour</td> <td></td> </tr> <tr> <td>Point 2:</td> <td></td> <td></td> <td></td> <td></td> <td>Calibration</td> <td><input type="text" value=""/></td> <td>Day</td> <td><input type="text" value=""/></td> <td>Hour</td> <td></td> </tr> <tr> <td>MB Type Factor:</td> <td colspan="5"><input type="text" value=""/></td> <td>1-Point Calibration Point</td> <td><input type="text" value=""/></td> <td colspan="3"><input type="checkbox"/> with Conc-0</td> </tr> </table>						Calibration Type:	<input type="text" value="6AB"/>	Formula:	<input type="text" value="Spline"/>	Counts:	<input type="text" value="2"/>						<Calibrator Parameters>						Calibrator	OD	Conc	Factor/OD range	Low	High	Slope Check <input type="text" value="None"/>					Point 1:	Saline			-0.1	2.5	Allowance Range Check <input type="checkbox"/> Reagent Blank <input type="checkbox"/> Calibration Advanced Calibration Operation <input type="text" value=""/> Interval (RB/ACAL) <input type="text" value=""/>					Point 2:	TruCal ASLO 1		*	-0.1	2.5	Point 3:	TruCal ASLO 2		*	-0.1	2.5	Point 4:	TruCal ASLO 3		*	-0.1	2.5	Point 5:	TruCal ASLO 4		*	-0.1	2.5	Point 6:	TruCal ASLO 5		*	-0.1	2.5						Point 7:											Point 8:											Point 9:											Point 10:											<Point Cal. For Master Curve>						Calibrator	OD	Conc	Low	High	Stability						Point 1:					Reagent Blank	<input type="text" value=""/>	Day	<input type="text" value=""/>	Hour		Point 2:					Calibration	<input type="text" value=""/>	Day	<input type="text" value=""/>	Hour		MB Type Factor:	<input type="text" value=""/>					1-Point Calibration Point	<input type="text" value=""/>	<input type="checkbox"/> with Conc-0																																
Calibration Type:	<input type="text" value="6AB"/>	Formula:	<input type="text" value="Spline"/>	Counts:	<input type="text" value="2"/>																																																																																																																																																																																																						
<Calibrator Parameters>																																																																																																																																																																																																											
Calibrator	OD	Conc	Factor/OD range	Low	High	Slope Check <input type="text" value="None"/>																																																																																																																																																																																																					
Point 1:	Saline			-0.1	2.5	Allowance Range Check <input type="checkbox"/> Reagent Blank <input type="checkbox"/> Calibration Advanced Calibration Operation <input type="text" value=""/> Interval (RB/ACAL) <input type="text" value=""/>																																																																																																																																																																																																					
Point 2:	TruCal ASLO 1		*	-0.1	2.5																																																																																																																																																																																																						
Point 3:	TruCal ASLO 2		*	-0.1	2.5																																																																																																																																																																																																						
Point 4:	TruCal ASLO 3		*	-0.1	2.5																																																																																																																																																																																																						
Point 5:	TruCal ASLO 4		*	-0.1	2.5																																																																																																																																																																																																						
Point 6:	TruCal ASLO 5		*	-0.1	2.5																																																																																																																																																																																																						
Point 7:																																																																																																																																																																																																											
Point 8:																																																																																																																																																																																																											
Point 9:																																																																																																																																																																																																											
Point 10:																																																																																																																																																																																																											
<Point Cal. For Master Curve>																																																																																																																																																																																																											
Calibrator	OD	Conc	Low	High	Stability																																																																																																																																																																																																						
Point 1:					Reagent Blank	<input type="text" value=""/>	Day	<input type="text" value=""/>	Hour																																																																																																																																																																																																		
Point 2:					Calibration	<input type="text" value=""/>	Day	<input type="text" value=""/>	Hour																																																																																																																																																																																																		
MB Type Factor:	<input type="text" value=""/>					1-Point Calibration Point	<input type="text" value=""/>	<input type="checkbox"/> with Conc-0																																																																																																																																																																																																			

Диапазон нормальных значений указан в соответствии с рекомендациями производителя реагентов. При использовании единиц измерения, отличающихся от приведенных, убедитесь, что значения стандартов, контрольных материалов, диапазонов нормальных значений и линейности метода введены в этих же единицах.

*-вводится из паспорта к калибратору (TruCal ASO, 5 уровней).

Контроль по TruLab Protein уровень 1 и уровень 2.