

# С -Реактивный белок, ЦРБ (CRP-FS)

ImmunoTD (Иммунотурбидиметрический)



Параметры для ввода в программу анализатора Architect c8000 / c4000

Configure assay parameters...

☒ General ☐ Calibration ☐ SmartWash ☐ Results ☐ Interpretation

Assay:  Type: Photometric Version: 1  
Number:  Assay availability:  Date: 18.02.2020  
Run controls for onboard reagents by:  Time: 11:01:08  
Operator: ADMIN

☒ Reaction definition ☐ Reagent / Sample ☐ Validity checks

Reaction mode:

Primary Secondary Read times  
Wavelength:  /  Main:  --   
Last required read:   
Absorbance range:  --  Color correction:  --   
Sample blank type:  Blank:  --

Assay 13 of 54

Done Cancel

?

Configure assay parameters...

☒ General ☐ Calibration ☐ SmartWash ☐ Results ☐ Interpretation

Assay:  Type: Photometric Version: 1  
Number:  Assay availability:  Date: 20.02.2020  
Run controls for onboard reagents by:  Time: 17:14:07  
Operator: FSE

☐ Reaction definition ☒ Reagent / Sample ☐ Validity checks

Reagent:   
Diluent:   
Diluent dispense mode:

R1 R2  
Reagent volume:    
Water volume:    
Dispense mode:

Dilution name	Sample	Diluted sample	Diluent	Water	Dilution factor	Default dilution
Standart	9.6				1:1.00	<input checked="" type="radio"/>
1:3	3.2				1:3.00	<input type="radio"/>
2x	18.0				2:1.00	<input type="radio"/>

Assay of

Done Cancel

?

# С -Реактивный белок, ЦРБ (CRP-FS)

ImmunoTD (Иммунотурбидиметрический)



Параметры для ввода в программу анализатора Architect c8000 / c4000

Configure assay parameters...

General Calibration SmartWash Results Interpretation

Assay: CRP Assay number: \*\*\* Date: 20.02.2020  
Calibration method: Spline Time: 17:14  
Operator: FSE

Calibrators Volumes Intervals Validity checks

Calibrator set: TruCalCRP

Calibrator level: Blank: Water Concentration: 0.0000

Cal 1: TruCalCRP1  
Cal 2: TruCalCRP2  
Cal 3: TruCalCRP3  
Cal 4: TruCalCRP4  
Cal 5: TruCalCRP5  
Cal 6: None

Replicates: 2 [Range 1 - 3]

Assay of

Done Cancel

?

Configure assay parameters...

General Calibration SmartWash Results Interpretation

Assay: CRP Assay number: \*\*\* Date: 20.02.2020  
Calibration method: Spline Time: 17:14  
Operator: FSE

Calibrators Volumes Intervals Validity checks

Calibrator: TruCalCRP

Calibrator level	Sample	Diluted sample	Diluent	Water
Blank: Water	9.6			
Cal 1: TruCalCRP1	9.6			
Cal 2: TruCalCRP2	9.6			
Cal 3: TruCalCRP3	9.6			
Cal 4: TruCalCRP4	9.6			
Cal 5: TruCalCRP5	9.6			
Cal 6: None				

Assay of

Done Cancel

?

# С -Реактивный белок, ЦРБ (CRP-FS)

*ImmunoTD (Иммунотурбидиметрический)*



Параметры для ввода в программу анализатора Architect c8000 / c4000

Configure assay parameters...

General Calibration SmartWash Results Interpretation

Assay: CRP Assay number: \*\*\* Date: 20.02.2020  
Calibration method: Spline Time: 17:14  
Operator: FSE

Calibrators Volumes Intervals Validity checks

Calibration intervals:

Full interval: 0 (hours)  
Adjust interval: 0 (hours)

Calibration type:

Adjust type: 2 Point  
Adjust level: Cal 1  
Default ordering type: Full

Assay of

Done  
Cancel

?

Единицы измерения - мг/л; знаков после запятой 2 (Configure result units).

Диапазон метода: 1-350 мг/л (Configure assay parameters\ Results).

Рефереренсный интервал: 0-5 мг/л (Configure assay parameters\ Results).  
При наличии ЛИС, не обязательно.

Диапазон нормальных значений указан в соответствии с рекомендациями производителя реагентов.

Калибратор TruCal CRP, набор из 5 уровней, нулевой - вода (или физраствор).

Контроль по TruLab protein уровень 1 и уровень 2.