

TruLab L Level 1

LOT 37094

Cat. No. 5 9020 99 10 065 (1 x 3 mL)

ВНИМАНИЕ!
 Уточняйте паспортные значения
 на www.diakonlab.ru или
 по тел.: 8-800-200-63-39



2027-09-30


 Updated sheets: www.diasys-diagnostics.com/service-area/value-sheets/

Value sheet version: 1

Short Name / Analyte	Method	Reagent Cat. No.	Value	Range	Unit
APOA1 Apolipoprotein A1 Аполипопротеин А1	Immunoturbidimetric Иммунотурбидиметрический	1 7102	118 1,18	94,3 - 141 0,943 - 1,41	mg/dL g/L
APOB Apolipoprotein B Аполипопротеин В	Immunoturbidimetric Иммунотурбидиметрический	1 7112	87,4 0,874	69,9 - 105 0,699 - 1,05	mg/dL g/L
CHOL Cholesterol Общий холестерин	CHOD-PAP	1 1300	156 4,02 1,56	134 - 177 3,46 - 4,59 1,34 - 1,77	mg/dL mmol/L g/L
HDLC HDL-Cholesterol ЛПВП Холестерин	Direct, colorimetric Прямой, колориметрический	1 3561	38,6 0,998 0,386	30,9 - 46,3 0,798 - 1,20 0,309 - 0,463	mg/dL mmol/L g/L
	CHOD-PAP Precipitant С осаждением	1 1300 + 1 3540	32,60 0,843 0,326	26,08 - 39,12 0,674 - 1,012 0,261 - 0,391	mg/dL mmol/L g/L
LDLC LDL-Cholesterol ЛПНП Холестерин	Direct, colorimetric Прямой, колориметрический	1 4131	100 2,59 1,00	80,3 - 120 2,08 - 3,11 0,803 - 1,20	mg/dL mmol/L g/L
NEFA Non-esterified fatty acids Свободные жирные кислоты	Enzymatic Ферментативный	1 5781	16,4 0,580 0,164	13,1 - 19,7 0,464 - 0,696 0,131 - 0,197	mg/dL mmol/L g/L
PL Phospholipids Фосфолипиды	Enzymatic Ферментативный			not assayed	
TRIG Triglycerides Триглицериды	GPO-PAP	1 5710	111 1,25 1,11	88,8 - 133 1,00 - 1,50 0,888 - 1,33	mg/dL mmol/L g/L



DiaSys Diagnostic Systems GmbH
 Alte Strasse 9 65558 Holzheim Germany
www.diasys-diagnostics.com