

TruLab L Level 2

LOT 37095

Cat. No. 5 9030 99 10 065 (1 x 3 mL)

ВНИМАНИЕ!
Уточняйте паспортные значения
на www.diaconlab.ru или
по тел.: 8-800-200-63-39



2027-09-30



Updated sheets: www.diasys-diagnostics.com/service-area/value-sheets/

Value sheet version: 1

Short Name / Analyte	Method	Reagent Cat. No.	Value	Range	Unit
APOA1 Apolipoprotein A1 Аполипопротеин А1	Immunoturbidimetric Иммунотурбидиметрический	1 7102	175 1,75	140 - 210 1,40 - 2,10	mg/dL g/L
APOB Apolipoprotein B Аполипопротеин В	Immunoturbidimetric Иммунотурбидиметрический	1 7112	137 1,37	109 - 164 1,09 - 1,64	mg/dL g/L
CHOL Cholesterol Общий холестерин	CHOD-PAP	1 1300	249 6,44 2,49	214 - 284 5,54 - 7,34 2,14 - 2,84	mg/dL mmol/L g/L
HDLC HDL-Cholesterol ЛПВП Холестерин	Direct, colorimetric Прямой, колориметрический	1 3561	57,9 1,50 0,579	46,3 - 69,5 1,20 - 1,80 0,463 - 0,695	mg/dL mmol/L g/L
	CHOD-PAP Precipitant С осаждением	1 1300 + 1 3540	55,07 1,424 0,551	44,06 - 66,08 1,139 - 1,709 0,441 - 0,661	mg/dL mmol/L g/L
LDLC LDL-Cholesterol ЛПНП Холестерин	Direct, colorimetric Прямой, колориметрический	1 4131	154 3,99 1,54	123 - 185 3,19 - 4,79 1,23 - 1,85	mg/dL mmol/L g/L
NEFA Non-esterified fatty acids Свободные жирные кислоты	Enzymatic Ферментативный	1 5781	51,1 1,81 0,511	40,9 - 61,3 1,45 - 2,17 0,409 - 0,613	mg/dL mmol/L g/L
PL Phospholipids Фосфолипиды	Enzymatic Ферментативный	1 5741	not assayed		
TRIG Triglycerides Триглицериды	GPO-PAP	1 5710	261 2,94 2,61	209 - 314 2,35 - 3,53 2,09 - 3,14	mg/dL mmol/L g/L



DiaSys Diagnostic Systems GmbH
Alte Strasse 9 65558 Holzheim Germany
www.diasys-diagnostics.com