

TruLab Urine

LOT 30530



Level 1

Cat. No. 5 9170 99 10 061 (6 x 5 mL)

ВНИМАНИЕ!
Уточняйте паспортные значения
на www.diaconlab.ru или
по тел.: 8-800-200-63-39



2022-09-30



Updated sheets: www.diasys-diagnostics.com/service-area/value-sheets/

Value sheet version: 1

Short Name / Analyte	Method	Reagent ID	Value	Range	Unit
ALB Microalbumin Микроальбумин	Immunturbidimetric Иммунотурбидиметрия	1 0242	7,66 0,766	5,52 - 9,81 0,552 - 0,981	mg/L mg/dL
AMY Amylase α-Амилаза	EPS G-7	1 0501	40,8 0,680	32,6 - 49,0 0,544 - 0,816	U/L μkat/L
BUN Urea N BUN	Urease UV Уреаза UV	1 3101	385 137	285 - 485 102 - 137	mg/dL mmol/L
CA Calcium Кальций	Arsenazo Арсеназо	1 1130	1,86 7,45	1,54 - 2,17 6,18 - 8,72	mmol/L mg/dL
	Phosphonazo Фосфаназо	1 1181	2,06 8,28	1,71 - 2,41 6,87 - 9,68	mmol/L mg/dL
CL Chloride Хлориды	BM 6010/C Indirect ISE	980132	55,0 195	47,3 - 62,7 168 - 222	mmol/L mg/dL
CREA Creatinine Креатинин	Jaffé, without compensation Яффе без компенсации	1 1711	69,4 6,14	52,8 - 86,1 4,67 - 7,61	mg/dL mmol/L
	Enzymatic PAP Ферментативный	1 1759	70,2 6,20	53,3 - 87,0 4,72 - 7,69	mg/dL mmol/L
D Specific Gravity Удельный вес	Refractometer Рефрактометр		1,011	1,000 - 1,020	
GLUC Glucose Глюкоза	GOD-PAP Глюкозооксидазный	1 2500 + 1 2550	24,9 1,38	19,4 - 30,4 1,08 - 1,69	mg/dL mmol/L
	HK/G6P-DH Гексокиназный	1 2511	26,0 1,44	20,3 - 31,7 1,13 - 1,76	mg/dL mmol/L
hCG Human Chorionic Gonadotropin ХГЧ (качественный)				Negative Отрицательный	
K Potassium Калий	BM 6010/C Indirect ISE	980131	30,7 120	25,5 - 35,9 99,7 - 141	mmol/L mg/dL
MG Magnesium Магний	Xylidyl blue Ксилидиновый синий	1 4610	2,29 5,58	1,70 - 2,89 4,13 - 7,03	mmol/L mg/dL
Na Sodium Натрий	BM 6010/C Indirect ISE	980130	74,2 171	64,6 - 83,8 148 - 193	mmol/L mg/dL
OSMO Osmolality Осмолярность	Freezing point Точка замерзания		398	319 - 478	mOsm/kg
pH	pH Meter pH-метр		6,21	5,50 - 6,90	
PO3 Phosphorus Фосфор	Molybdate UV, substrate start Молибдат UV, запуск субстратом	1 5211	25,9 8,37	20,2 - 31,6 6,53 - 10,2	mg/dL mmol/L
TP Total Protein Общий белок	Pyrogallol red Пирогаллоловый красный	1 0210	7,50 75,0	5,70 - 9,30 57,0 - 93,0	mg/dL mg/L
UREA Urea Мочевина	Urease UV Уреаза UV	1 3101	824 137	610 - 1038 102 - 173	mg/dL mmol/L
UA Uric Acid Мочевая кислота	ТВНВА TOOS	1 3021 1 3001	7,90 470 8,48 504	5,85 - 9,95 348 - 592 6,28 - 10,7 373 - 636	mg/dL μmol/L mg/dL μmol/L



DiaSys Diagnostic Systems GmbH
Alte Strasse 9 65558 Holzheim Germany
www.diasys-diagnostics.com