

Всем пользователям анализаторов биохимических
автоматических CA-270 и CA-400

Мы, компания Фуруно Электрик Ко., Лтд., настоящим заявляем, что для работы на биохимических анализаторах CA-270 и CA-400 требуется специальная подготовленная вода.

Согласно техническим характеристикам FURUNO CA-270 потребляет не более 5 литров в час, в свою очередь CA-400 потребляет не более 18 литров в час.



Mitsuhiro Takeda

General Manager, ME Marketing Department, ME Business Unit
Furuno Electric Co., Ltd.

Date: 06 October, 2020

To all users of clinical chemistry analyzers
of CA-270 and CA-400

We, Furuno Electric Co., Ltd., herewith declare that the biochemical analyzers CA-270 and CA-400 require the special prepared water to work. According to the technical specifications, the CA-270 analyzer consumes no more than 5 liters per hour, while the CA-400 analyzer consumes no more than 18 liters per hour.



Mitsuhiro Takeda
General Manager, ME Marketing Department, ME Business Unit
Furuno Electric Co., Ltd.
Date: 06 October, 2020

